

# Short time exposure (STE) test protocol

## Index

### Description of the test methods

1. Cells
2. Preculture
3. Selection of solvents for test substance preparation
4. Preparation of test substance solutions
5. Experimental procedures

### Data analysis

6. Calculation of cell viability
7. Criteria for irritant and non-irritant classifications
8. Rank classification of eye irritation potential in STE test


### Additional information

9. Culture media and materials

## Test methods

Outline of the STE test procedures is described in Fig. 1.

Fig. 1 STE test methods



MTT: Methyl Thiazolyl Diphenyl-Tetrazolium Bromide

### 1. Cells

- Use SIRC cells purchased from ATCC (catalog No. CCL60, lot No. 3981569). The cells should be used between 3 weeks and 3 months after the start of cultivation or within 25 passages.
- Culture SIRC cells at 37°C under 5% CO<sub>2</sub> and humidifying in a culture flask containing Eagle's MEM supplemented with 10% (v/v) fetal bovine serum, 2 mM L-glutamine, 50-100 unit/mL penicillin and 50-100  $\mu$ g/mL streptomycin. Dissociate the cells that have become sub-confluent in the culture flask using trypsin-EDTA solution and subculture in a culture flask or inoculate into 96-well plates.


## 2. Preculture

- Prepare cell suspensions of  $3.0 \times 10^4$  cells/mL or  $1.5 \times 10^4$  cells/mL in culture medium. Preculture the cells for 4 days at 37°C under 5% CO<sub>2</sub> and humidifying if 200 µL each of the  $3.0 \times 10^4$  cells/mL suspension is to be added to each well of the 96-well plate; for 5 days if the  $1.5 \times 10^4$  cells/mL suspension is to be used for this procedure.

## 3. Selection of solvents for test substance preparation

A flowchart for the procedures from selection of solvents to preparation of substance is shown in Fig. 2.

**Fig. 2 Solvent selection and substance preparation**



- First, prepare 5% (w/w) solution of a test substance using saline as solvent and observe the dissolution pattern of the substance. If the substance dissolves or homogeneously disperses <sup>Note 1), 2)</sup>, choose saline as the solvent for this substance.
- When the substance does not dissolve or homogeneously disperse in saline, try saline containing 5% (w/w) DMSO. If the substance dissolves or homogeneously disperses, choose saline containing 5% DMSO as the solvent for this substance.

- When the substance does not dissolve or homogeneously disperse in saline containing 5% (w/w) DMSO, prepare this substance using mineral oil as solvent. Prepare 5% (w/w) solution of the substance using mineral oil as solvent and observe the dissolution pattern of the substance. If the substance dissolves or homogeneously disperses, choose mineral oil as the solvent for this substance. If the substance does not dissolve or homogeneously disperse in mineral oil, do not evaluate this substance in this test.

Note 1): A condition in which a substance is dispersed homogeneously and remains dispersed for at least 5 minutes

Note 2): Vortex, sonication or proper warm up can increase the solubility.

#### 4. Preparation of test substance solutions

For assessment of the transferability, test substance is evaluated with 3 concentrations of 5%, 0.5% and 0.05%. For the main test in which blinded test substance is evaluated, 3 concentrations of test substance of 5%, 0.5% and 0.05% are prepared but two concentrations (5% and 0.05%) are tested.

- Prepare sample solutions of the test substances using the solvents selected as described in Section 3.
- Weigh the test substance in a screw tube and prepare 5% (w/w) solution using the selected solvent. Then, prepare 0.5% (v/v) and 0.05% (v/v) solutions by serial 10-fold dilutions of the original 5% (w/w) solution.
- Use SDS 0.01% (v/v) saline solution as the positive control. Prepare SDS 1% (w/w) saline solution, and then prepare SDS 0.1% (v/v) and SDS 0.01% (v/v) saline solutions by serial 10-fold dilutions starting from the original SDS 1% (w/w) solution.
- Use the solvents used to prepare sample solutions as the solvent controls.

#### 5. Experimental procedures

- Prior to initiating treatment conditions, familiarize yourself with the treatment diagram described in Figure 3.
- Preculture the cells to confluency in a 96-well plate.
- Take 0.6 mL (for 0.2 mL  $\times$  3 wells) of the prepared test sample using a 1 mL disposable syringe.


- Tilt the plate and remove culture medium from each well by suction using a suction tube (This can be done more easily by connecting a Pasteur pipette to tip of the suction tube. Be careful not to touch the well bottom with tip of the Pasteur pipette) .
- Add test samples into wells. Start stopwatch when the first sample is added and then add test samples into their designated wells at a rate of 3 wells every 7-10 seconds (see Fig. 3) and record the elapsed time by this.
- Exactly five minutes later, sequentially remove the test samples from all wells by suction at a rate of 3 wells every 7-10 seconds. Stop stopwatch and record the total working time by this.
- Carefully add 200  $\mu$ L each of PBS into the wells from which the test samples have been removed, and then remove it by suction. Repeat this manipulation twice to wash the inner surfaces of the wells (use of 8-channel pipettor, etc. will facilitate this wash).
- Thoroughly remove PBS by suction so that not a trace of it will remain in any of the wells.
- Add 200  $\mu$ L of 0.5 mg/mL MTT (CAS No. 298-93-1) dissolved in culture medium <sup>Note 3)</sup> to each well and incubate in an incubator (37°C, 5% CO<sub>2</sub>) for 2 hours.

Note 3): In the case of 21 mg MTT, add 42 mL of MEM culture medium, shake well and let the mixture stand for 20 minutes. Shake the mixture well, place it in a 20 mL syringe and filtrate through a 0.45  $\mu$ m microfilter.

Protect MTT from light after weighing until immediately before use.


**Fig. 3 Exposure and Evaluation of Treated Cultures**

**Evaluation of standard substances**




Diagonal lines: PBS to prevent drying (200  $\mu$ L/well)  
 1-3'': Test samples (3 compounds, 3 doses)  
 P: Positive control (0.01% SDS)  
 S: Saline  
 C: Procedural control (culture medium with cells)  
 B: Blank (culture medium without cells)

1. Remove culture medium.
2. Add test samples starting at the wells marked 1. Start stopwatch and then add test sample solutions into designated wells every 7-10 seconds and record the elapsed time by this.
3. Five minutes later, suction off the test sample solutions every 7-10 seconds starting at the wells marked 1. Stop stopwatch and record the total working time by this.
4. Wash the wells and add MTT.




**Evaluation of test substances**



Diagonal lines: PBS to prevent drying (200  $\mu$ L/well)  
 1-2'': Test samples (3 compounds, 2 doses)  
 P: Positive control (0.01% SDS)  
 S: Saline  
 D: Saline containing 5% DMSO  
 M: Mineral oil  
 C: Procedural control (culture medium with cells)  
 B: Blank (culture medium without cells)

1. Remove culture medium from all wells (except wells for blank).
2. Add test samples starting at the wells marked 1. Start stopwatch and then add test sample solutions into designated wells at a rate of 3 wells every 7-10 seconds and record the elapsed time by this.
3. Five minutes later, sequentially suction off the test sample solutions at a rate of 3 wells every 7-10 seconds starting at the wells marked 1. Stop stopwatch and record the total working time by this.
4. Wash the wells and add MTT.



- Take out the plate from incubator and discard the liquid in the wells. Remove the remaining drops completely by lightly tapping the plate onto a paper towel.
- Add 200  $\mu$ L of 0.04N HCl-isopropanol (isopropanol 96 mL + 1N HCl 4 mL) to each well and let the plate stand for 60 minutes.
- Make the extracts in the wells uniform, set the plate onto a plate reader and measure absorbance at 570 nm wavelength.

#### Data analysis

##### 6. Calculation of cell viability

- For each test sample, calculate cell viability (%) as the ratio of the optical density of test sample to the optical density of solvent control.

$$\text{Cell viability} = \frac{\text{optical density of test sample (optical density of test sample - optical density of blank)}}{\text{optical density of solvent control (optical density of solvent control - optical density of blank)}} \times 100$$

- In each assay, use 3 wells for each sample concentration and use the average of 3 optical density measurements to calculate cell viability.
- The cell viability of test sample concentration is treated as “zero”, when the average of 3 well optical densities is less than optical density of solvent control and the cell viability shows minus.
- Perform 3 independent assays\*1 for each sample concentration and use the average of 3 assays as the final cell viability.

#### Remark

\*1: Each test sample concentration is evaluated with 3 wells on a plate. In case the following 1) to 3) criteria are satisfied on the plate, the average cell viability of 3 wells is treated as an acceptable assay independently.

#### Criteria for acceptable test

- 1) Optical density of the culture medium procedural control (C in Fig. 3) at least 0.3 after subtraction of blank optical density
- 2) Cell viability for the solvent control at least 80% relative to that for the culture medium procedural control
- 3) Cell viability for the positive control (0.01% SDS) within a range of 21.1 to 62.3 % (within a range of mean cell viability  $41.7 \% \pm 10.3 \times 2$  standard deviations)

4) Standard deviation of the final cell viability derived from 3 independent assays less than 15 % for all test sample concentrations\*2

Remark

\*2: In case the standard deviation of cell viability is more than 15 %, additionally other 3 assays are operated and the above criteria 4) is evaluated with independent results for additional 3 assays.

7. Criteria for irritant and non-irritant classifications

- In the STE test, cell viability is the designated endpoint after a 5-minute exposure to test sample. Substances demonstrating higher than 70% cell viability are classified as non-irritant (NI) and those demonstrating a cell viability equal to or less than 70% are classified as irritant (I).

Classification	Cell viability
Non-irritant	> 70%
Irritant	≤ 70%

8. Rank classification of eye irritation potential in STE test

- The procedure for rank classification of eye irritation potential in the STE test is shown in Fig. 4.
- With the 5% test sample solution, a score of 0 is given when cell viability is higher than 70% and score of 1 is given when cell viability is not higher than 70%. For 0.05% test sample solution, score of 1 is given when cell viability is higher than 70% and a score of 2 is given when cell viability is not higher than 70%. Next, the scores are added up for the 5% and 0.05% test sample solutions and final value is compared to the rank scale. The substance is classified according the following scale: Rank 1 = score of 1 = minimally irritant; rank 2 = score of 2 = moderately irritant; rank 3 = score of 3 = severely irritant.

**Fig. 4 Rank classification of eye irritation potential in STE test**

**STE irritation score**

Test sample solution 5%	Score	Test sample solution 0.05%	Score
If CV > 70%,	0	If CV > 70%,	1
If CV ≤ 70%,	1	If CV ≤ 70%,	2

CV: Cell viability



**5% score + 0.05% score = STE rank**

STE rank	Eye irritation potential
1	Minimally irritant
2	Moderately irritant
3	Severely irritant

## Additional Information

### 9. Culture media and materials

#### 9-1. Preparation of culture medium

- Add the following materials to 500 mL of E-MEM:

-FBS 56 mL

-L-glutamine 6 mL

-Penicillin-streptomycin 6 mL

#### 9-2. Precautions in preparing culture medium

- If white turbidity is seen in L-glutamine solution, dissolve it in a constant-temperature bath.
- Prior to the preparation of penicillin-streptomycin solution, lightly shake the container bottle to mix the solution.
- Use disposable pipettes only once throughout the preparation procedure (use two 25 mL pipettes for FBS).

#### 9-3. STE test-related agent and materials

Item	Maker	Item code	Note
SIRC cells*	ATCC	CCL-60	Lot.3981569
Eagle's MEM medium*	Sigma-Aldrich	M2279-500ML	
Fetal bovine serum*	GIBCO	not yet decided	Lot. not yet decided
L-glutamine 200 mM*	GIBCO	25030-081	
Penicillin-streptomycin*	GIBCO	15140-122	0.85% NaCl Penicillin 10,000U Streptomycin 10mg
MTT*	Sigma-Aldrich	M-2128	
Saline for injection*	Otsuka Pharmaceutical	081517	Japanese Pharmacopoeia, 20mL, 50 tubes
Dimethyl sulfoxide (DMSO)*	SIGMA-Aldrich	154938	
Mineral oil*	SIGMA	M5310	Mineral oil, Embryo tested; sterile-filtered
PBS tablet	Takara Bio	T900	
Trypsin-EDTA	Sigma-Aldrich	T-3924	Porcine trypsin 0.5 g EDTA-4 Na 0.2 g in HBSS 1 L

2-propanol	Kanto Chemical	32435-01	
1 mol/L hydrochloric acid (1N)	Kanto Chemical	18591-08	
Sumilon cell culture flask (middle size)	Sumilon	MS-21250	Tissue Culture Treated Polystyrene Sterile
96-well microplate (flat bottom)	Corning	3997	Tissue Culture Treated Polystyrene Sterile
Tissue culture dish	Corning	430165	
Sumilon serum tube	Sumitomo Bakelite	MS-4601X	
Cell scraper	Sumilon	MS-93170	
Terumo syringe	Terumo	SS-20ESZ	
Microfilter Millex-HA filter unit	Millipore	SLHA033SS	MF-Millipore (cellulose mixed ester)
Falcon tube	Becton Dickinson (Falcon)	2070	
15 mL tube	Asahi Techno Glass	2314-015	
Sumilon pipette	Sumilon	MS-66050	
Disposable pipette	Sumilon	MS-66250	
Eppendorf combi-tip plus (Sterilized)	Eppendorf	30069455	
SCC screw tube bottle (white) No. 2	As One	7-2110-04	

\*Use only product of the same specification. For other reagents and materials, equivalents regardless of makers thereof are also acceptable.