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## 1. Purpose

This test method is used to measure cytotoxicity of chemicals using SIRC cell and to discriminate between non ocular irritant (NI) and ocular irritant (Others) as a bottom up approach<sup>1-3)</sup>. The *in vivo* standard for the assessment is based on the classification of the Globally Harmonized System of Classification.

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## 2. The principle of SIRC cytotoxicity test

The cytotoxicity has been known as an index for evaluating the eye irritation by chemicals, as the corneal damage that has a greater impact on the total eye irritation is related to that of the corneal epithelium cell<sup>4)</sup>. The cytotoxicity test is reported to be useful for identifying the non ocular irritants that has almost no effect on the cornea. The cell line used comes from rabbit cornea. The application time of 72hr is set with consideration of time (within 72hr) from ocular instillation to maximal eye irritation for general chemicals except acid or alkali, on the basis of *in vivo* data at the previous research project<sup>5)</sup>.

The SIRC cytotoxicity test is based on the measurement of viable cells stained by crystal violet, which penetrates via cell membrane and stains biological macromolecules. The crystal violet staining method can be used for many cultured cells and can produce the relatively invariable results<sup>5-9)</sup>. The relative control is additionally used to obtain the invariable results<sup>10,11)</sup>. Moreover, the operation is simple and easy, and the tested microplate can be stored. The results can be confirmed by using the stored microplates in any time. No other method can match it.

On the other hand, the disadvantage of this method is to be confined to test substances which are solved or uniformly suspended in the medium.

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## 3. Materials

### 3.1. Cell

The SIRC cell (Statens Seruminstitut rabbit corneal cell: ATCC NO. CCL-60), the cell line from the rabbit cornea is obtained from ATCC (American Type Culture Collection). The cells can be frozen and kept in liquid nitrogen. The cells should be confirmed the absence of mycoplasma (e.g., Venor GeM Mycoplasma Detection Kit, Minerva Biolabs GmbH, 11-1025). The cells should be used during 3 months after the start of cultivation (e.g. When the cell culture starts at passage number of 435 and is maintained by passage every four days, it

should be used within passage number of 470). They should be checked on the basis of the quality control of 4.7.

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Minimum Essential Medium (MEM)

Fetal Bovine Serum (FBS)

・Inactivated fetal bovine serum should be used. Inactivation is performed by 56°C, 30 min in the water bath. After fall in temperature, the serum is taken 56 mL or 28 mL in each tube. The serum is stored at -70°C or less or at -20°C or less.

Penicillin/Streptomycin/Amphotericin B (P/S/F) solution

(Antibiotic-Antimycotic 100x ; GIBCO BRL)

Phosphate-Buffered Saline (-) (PBS(-))

0.25% (w/v) Trypsin (1mmol/L EDTA·4Na)

Dimethyl Sulfoxide (DMSO; CAS Number 67-68-5)

・the weighing can be converted to volume on the basis of specific gravity.

Ethanol (CAS Number 64-17-5)

・the weighing can be converted to volume on the basis of specific gravity.

Crystal Violet (CAS Number 548-62-9)

Methanol (CAS Number 67-56-1)

Sodium Dodecyl Sulfate (SDS; CAS Number 151-21-3)

Triethanolamine (CAS Number 102-71-6)

・It should be used that of Purity >= 98.0%.

Hydrochloric Acid (CAS Number 7647-01-0)

Sodium Hydroxide (CAS Number 1310-73-2)

The example of the manufacturer and so on of the reagent is shown in the Table 1.

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### 3.5. Medium

MEM is supplemented with 10% FBS (inactivated) and 1% appropriate antibiotic (P/S/F solution). For example, 500 mL of MEM is supplemented with 56 mL of FBS and 5.6 mL P/S/F. At this time, the concentrations of the antibiotics are Penicillin 100 U/mL, Streptomycin 100 µg/mL and Amphotericin B 250 ng/mL, respectively.

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### 3.6. Crystal violet solution

0.4% Crystal violet solution is prepared using methanol.

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

### 3.7.1. Examination of stability for the substance in the medium

The solubility of each substance in the medium should be confirmed in advance, using the procedure shown in Fig. 1. The solubility or suspensibility of the substance is examined at the concentration of 10,000 µg/mL (1 w/v%) in the medium. It is solved or uniformly suspended using vortex mixer, waterbath and sonicator when it finds necessary. When it is not suspended uniformly, the solubility to DMSO and ethanol is examined as a next step. The substance and the selected solvent are mixed at the concentration of 10,000 µg/mL (1 w/v%) in the medium, respectively. The suspensibility is examined after sonication and so on. When it is not suspended uniformly at the concentration of 10,000 µg/mL (1 w/v%), the suspensibility is examined using the substance of 5,000 µg/mL (0.5 w/v%) containing the selected solvent of 10,000 µg/mL (1 w/v%) in the medium. The substance that is not suspended homogeneously is judged as an inapplicable substance for this test, and has no test. These judges are performed macroscopically.

Figure 1. Flow chart of examination of stability for the substance in the medium

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### 3.7.2. Preparation of test substances

The test substance is solved or uniformly suspended at appropriate concentration by the

procedure described in 3.7.1. The final maximal concentrations of the substances and the recommended solvents are 5,000 µg/mL (0.5 w/v%), respectively, after diluted by the medium containing the SIRC cells as described in 4.3. The final maximal concentration of the substances and recommended solvents are 2,500 µg/mL (0.25 w/v%) and 5,000 µg/mL (0.5 w/v%), respectively, when the low concentration of the substance is selected. Furthermore, the data of the well with precipitation or so on at anytime after the mix of the substance and the cell, especially after 72 hr incubation, should be rejected for unsuitable suspension.

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### 3.7.3. Dilution of substance solution

The dilution step of substance should be eight by a common ratio of two. The test of the same concentration should be performed in duplicate on the plate.

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## 3.8. Reference substance

### 3.8.1. Positive control substance

SDS should be used as a positive control. The SDS solution is prepared at the concentration of 1,000 µg/mL in the medium.

### 3.8.2. Relative control substance

Triethanolamine should be used as a relative control. The triethanolamine solution is prepared at the concentration of 10,000 µg/mL in the medium.

### 3.8.3. Negative control substance

Medium, 10,000 µg/mL DMSO-medium solution or 10,000 µg/mL ethanol-medium solution

should be used as a negative control. It should be selected on the basis of what kind of solvent solve or suspend test substance.

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##### 4.1. Passage culture of cell

###### ・cell culture

- (1)The SIRC cell is cultured with MEM supplemented 10% FBS, 1% P/S/F and so on at 37 °C in a humidified incubator containing 5% (v/v) CO<sub>2</sub> in air. The concentrations of the penicillin, streptomycin and amphotericin B are 100 U/mL, 100 µg/mL and 250 ng/mL, respectively.
- (2)After the medium is removed from culture flask, the SIRC cells are rinsed 2 times by PBS(-)

10 mL for removing serum as a trypsin inhibitor.

- (3) After PBS(-) is removed, 0.25% trypsin solution (1.5-2 mL) is added contacting the all cells in the culture flask.
- (4) After 0.25% trypsin solution is removed halfway, the cells are incubated for 2-3 minutes at 37°C. The cells are detached from the surface of the flask, tapping it. They are collected by the appropriate amount of the MEM (10% FBS) from the culture flask. After counting the cells, the cell suspension is prepared at the density of  $6-8 \times 10^5$  cells in 15-30 mL of the medium. The passage culture is performed using it.

•Freeze preservation of cell

- (1) The medium contained 10% DMSO is prepared as a medium for the freeze preservation of cell. The commercial available freeze preservative solution (e.g., Juji Field Inc., Cell banker 1 or 2) can be used. Both serum type and non serum type can be used.
- (2) The cell of  $1 \times 10^6$  cells/mL is add to stock tube. It is froze slowly. For example, the freeze is performed at 5 minutes in the ice, 50 minutes at about -20°C and 12 hours in about -70°C and in the liquid nitrogen in order. The commercial available freezing vessel (e.g., Nihon Freezer Co. LTD., BICELL) can be used. In this case, the vessel containing the tube for the cell is frozen at about -70°C and in the liquid nitrogen after 3 hours - overnight.
- (3) The tube containing the cell is preserved in the liquid nitrogen.

•Melting of cell

- (1) The stock tube is dipped in the water of 37°C for melting.
- (2) The cell suspension added the medium of 10 mL is centrifuged at 1,000 rpm x5 min.
- (3) After removing the supernatant, the cell suspension is prepared by the medium.

•The freeze preservative cell is passaged one or more times and should be confirmed the appropriate growth.

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#### 4.2. Preparation of cell suspension

- (1) After the medium is removed from culture flask, the SIRC cells are rinsed 2 times by PBS(-) 10 mL for removing serum as a trypsin inhibitor.
- (2) After PBS(-) is removed, 0.25% trypsin solution (1.5-2 mL) is added contacting the all cells in the culture flask.
- (3) After 0.25% trypsin solution is removed halfway, the cells are incubated for 2-3 minutes at 37°C.
- (4) The cells are detached from the surface of the flask, tapping it.
- (5) The cells are collected by the appropriate amount of the MEM (10% FBS) and pepetting.
- (6) After counting the cells, the cell suspension is prepared at the density of  $2 \times 10^5$  cells/mL.

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#### 4.3. Application of test substance

- (1) PBS(-), negative control, and dilution series of test substance, positive control and relative control are prepared in the 96 well microplate. The layout of the microplate is shown in fig. 1.
- (2) One tenth mL of the  $2 \times 10^5$  cells/mL cell suspension is added to the wells as shown in fig.2.
- (3) Microplate sealer is used to avoid the effect of volatile toxicant. Moreover, wrap film can be used. The five measurements (4.7.(1)-(5)) of the quality control should be used for checking wether the volatile substance has an effect on other wells. The criterion of the toxic effect is the same as that of the quality control. When the volatile substance has an effect on other wells, the retest should be performed using dilution.
- (4) The microplate added the test substance and the cell suspension is still standing for 20 minutes in the clean bench. The cell is adhered to the bottom of the microplate. After that, the microplate is moved to the CO<sub>2</sub> incubator.
- (5) The microplate is incubated for about 72 hr at the condition of 37°C and 5% CO<sub>2</sub>,

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#### 4.4. Crystal violet staining

- (1)After the incubation, the medium containing the test substances is removed by gentle turnover of the microplate.
- (2)PBS(-) of 200  $\mu$ L is added. After gentle shake, PBS(-) is removed by the turnover of the microplate. This procedure is repeated twice.
- (3)One hundred  $\mu$ L of crystal violet methanol solution is added to each well of the microplate, and stain the cells for 30 minutes.
- (4)After the staining, the crystal violet methanol solution is removed by gentle turnover of the microplate. It is washed by tap water sufficiently, and is blotted water with the paper towel.
- (5)After drying sufficiently, the absorbance at 588 nm is measured with an automatic microplate reader. The measurement at the near wavelength that the equivalency was demonstrated can be used.

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#### 4.5. Calculation of IC<sub>50</sub>

The absorbance of control wells, which contained no test substance, is regarded as 100%, and the percentage absorbance for each well is calculated. The concentration at which the growth of cells was inhibited to 50% of the control (IC<sub>50</sub>) is obtained from the next formula using two concentrations around the predicted concentration of 50% cell viability.

$$\text{LogIC}_{50} = [(50-y_1)\text{log}x_2 - (50-y_2)\text{log}x_1] / (y_2 - y_1)$$

(x<sub>1</sub> : Low concentration, x<sub>2</sub> : High concentration, y<sub>1</sub>: cell viability at low concentration, y<sub>2</sub> : cell viability at high concentration, Log means common logarithm.)

If the cell viability is >50% at the maximal concentration of 5,000 μg/mL, the result of the test substance is IC<sub>50</sub>>5,000 μg/mL. Also, if the cell viability is <50% at the concentration of 39.1 μg/mL, the result of the test substance is IC<sub>50</sub><39.1 μg/mL.

If the multiple concentrations showing 50% cell viability were obtained from one substance, the lowest IC<sub>50</sub> should be adopted.

In a spreadsheet software (Excel), the cell viability value is rounded to the nearest tenths.

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#### 4.6. Evaluation

The eye irritation of the test substance is evaluated by using triethanolamine as a relative control (see annex 1). The classification of the eye irritation of triethanolamine (100%) is “Not Category” of GHS standard. For the GHS standard, the test substance is judged as negative (Not Category) when the IC<sub>50</sub> is higher than or equal to that of triethanolamine, and is judged as positive (Category 1 or 2) when the IC<sub>50</sub> is lower than that of triethanolamine. If the two decisions are different, the two results with the same decision are adopted after the third test. If the confirmation of the deviation is needed with three results, the test is repeated three times.

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#### 4.7. Quality control

The quality control of the SIRC cytotoxicity test is performed by the six measurements. It needs to satisfy the next criterion. If it is not satisfy the criterion, the retest of the microplate should be performed. Especially, if the toxic effect of the volatile substance causes unacceptable values, the retest should be performed at the lower concentration.

- (1)The absolute OD obtained from the negative control is index of normal proliferation of SIRC cell seeded at the concentration of  $1 \times 10^4$  cells/well and incubated 72hr. The mean OD of the negative control (the right and left wells) should be  $> 0.4$  for complete test.
- (2) Sodium dodecyl sulfate (SDS) is used as a positive control. The  $IC_{50}$  of SDS should be 77.7-258.7  $\mu\text{g}/\text{mL}$  when it is tested by the standard protocol. It needs for complete test.
- (3) Triethanolamine is used as a relative control (see appendix 1). The  $IC_{50}$  range of triethanolamine should be 1,000-2,500  $\mu\text{g}/\text{mL}$  when it is tested by the standard protocol (see appendix 2 ). It needs for complete test.
- (4) The difference between two dilution series of the substance should be confirmed. The  $IC_{50}$ s of the first series and the second series should be within  $\pm 20\%$  of the mean  $IC_{50}$  of two series(the mean of the first  $IC_{50}$  and the second  $IC_{50}$ ), respectively. It needs for complete test. If the  $IC_{50}$  is lower than 39.1  $\mu\text{g}/\text{mL}$ , it is calculated using the  $IC_{50} = 39.1 \mu\text{g}/\text{mL}$ . If the

IC<sub>50</sub> is higher than 5000μg/mL, it is calculated using the IC<sub>50</sub>=5000μg/mL. Treatment of IC<sub>50</sub> expressed with inequality sign is performed in the same manner.

- (5) The difference between left and right wells of the negative control should be confirmed to check the systematic quality. The mean ODs of left and right wells should be within  $\pm 15\%$  of the mean OD of all negative control wells (the mean  $\pm 15\%$ ), respectively. It needs for complete test.
- (6) The difference between the two test results adopted should be confirmed to check the equality of them. The IC<sub>50</sub> values of two tests of positive control (SDS) should be lower or equal to twice. (The high value/the low value  $\leq 2$ )

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- 1) Scott, L. et al., *Toxicol. In Vitro*, 24(1), 1-9 (2010).
- 2) Hagino, S. et al., *ATLA*, 36, 641-652 (2008)
- 3) Hagino, S. et al., *ATLA*, 38, 139-152 (2010)
- 4) Itagaki, H. et al., *Toxicol. in Vitro*, 5, 139-143 (1991).
- 5) Ohno, Y. et al., *Toxicol. in Vitro*, 13, 73 (1999).
- 6) Saotome, K. et al., *Toxicol. in Vitro*, 3, 317-321 (1989).
- 7) Itagaki, H., *AATEX*, 3, 182-190 (1995).
- 8) Ohno, Y et al., *AATEX*, 3, 123 (1995).
- 9) Tani, N., *Toxicol. in Vitro*, 13,175 (1999).
- 10) 「代替法を用いて化粧品原料の眼刺激性を評価するにあたっての指針(厚生科学研究班の作成した案)」*AATEX*,5,Suppl., Guideline Draft1-3 (1998).
- 11) Ohno, Y., *ATLA*, 32, Supplement 1, 643-655, 2004.

## 5.References

- 1) Scott, L. et al., *Toxicol. In Vitro*, 24(1), 1-9 (2010).
- 2) Hagino, S. et al., *ATLA*, 36, 641-652 (2008)
- 3) Hagino, S. et al., *ATLA*, 38, 139-152 (2010)
- 4) Itagaki, H. et al., *Toxicol. in Vitro*, 5, 139-143 (1991).
- 5) Ohno, Y. et al., *Toxicol. in Vitro*, 13, 73 (1999).
- 6) Saotome, K. et al., *Toxicol. in Vitro*, 3, 317-321 (1989).
- 7) Itagaki, H., *AATEX*, 3, 182-190 (1995).
- 8) Ohno, Y et al., *AATEX*, 3, 123 (1995).
- 9) Tani, N., *Toxicol. in Vitro*, 13,175 (1999).
- 10) Guidance for evaluation of eye irritation of cosmetic ingredients using alternative method (Draft document by the study team supported by Ministry of Health and Welfare), *AATEX*,5,Suppl., Guideline Draft1-3 (1998).
- 11) Ohno, Y., *ATLA*, 32, Supplement 1, 643-655, 2004.

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°C	Degrees Centigrade
ATCC	American Type Culture Collection
DMSO	Dimethyl Sulfoxide
EPA	United States Environmental Protection Agency
FBS	Fetal Bovine Serum
GHS	Globally Harmonized System of Classification and Labelling of Chemicals
IC <sub>50</sub>	50% Inhibitory Concentration
JaCVAM	Japanese Center for the Validation of Alternative Methods
MEM	Minimum Essential Medium
NI	Non Irritant
OD	Optical density
PBS(-)	Phosphate-Buffered Saline (-)
SDS	Sodium Dodecyl Sulfate
SIRC cell	Statens Seruminstitut Rabbit Corneal Cell
SIRC-CVS	Statens Seruminstitut Rabbit Cornea-Crystal Violet Staining

## 6.List of abbreviations and acronyms

°C	Degrees Centigrade
ATCC	American Type Culture Collection
DMSO	Dimethyl Sulfoxide
EPA	United States Environmental Protection Agency
FBS	Fetal Bovine Serum
GHS	Globally Harmonized System of Classification and Labelling of Chemicals
IC <sub>50</sub>	50% Inhibitory Concentration
JaCVAM	Japanese Center for the Validation of Alternative Methods
MEM	Minimum Essential Medium
NI	Non Irritant
OD	Optical density
PBS(-)	Phosphate-Buffered Saline (-)
SDS	Sodium Dodecyl Sulfate
SIRC cell	Statens Seruminstitut Rabbit Corneal Cell
SIRC-CVS	Statens Seruminstitut Rabbit Cornea-Crystal Violet Staining

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#### 7.About the revision of this protocol

(1) The revision of ver.1 - ver.1.71 that is the same as ver.1.71j and ver.1.71e, is shown by the green character in the ver 2.13. The protocol of the ver.1 was used for evaluating 68 chemicals at Shiseido Research Center in 2009-2010 and was subjected to the peer review by JaCVAM.

(2) The revision after ver.1.71 is shown by the blue character in the ver.2.13.

(3) The revision from ver.2.07 to ver.2.08:

In 4.7.(4), “The difference between two dilution series of the substance should be confirmed. The IC50s of the first series and the second series should be within + 20% of the mean IC50 of two series (the mean  $\pm$  20%), respectively.” was changed to “The difference between two dilution series of the substance should be confirmed. The IC50s of the first series and the second series should be within  $\pm$  20% of the mean IC50 of two series(the mean of the first IC50 and the second IC50), respectively.”

(4) The revision form ver.2.08 to ver.2.09:

In 4.3.(3), “The five measurements (4.7.(1)-(5)) of the quality control should be used for checking whether the volatile substance has an effect on other wells. The criterion of the toxic effect is the same as that of the quality control. When the volatile substance has an effect on other wells, the retest should be performed using dilution.”

(5) The revision from ver.2.09 to ver.2.11:

“SIRC cytotoxicity test” was changed to “SIRC-CVS cytotoxicity test”. “Table 1” of 3.7.1. was changed to “Figure 1”.

(6) The revision from ver.2.11 to ver.2.12:

The version was added to title. SIRC-CVS was added to list of abbreviations and acronyms of 6. “Figure 1. Layout of 96 well microplate” was changed to “Figure 2. Layout of 96 well microplate”. “Figure 2. Addition of cell suspension” was changed to “Figure 3. Addition of cell suspension”

(7)The version from ver.2.12 to ver.2.13:

In 3.7.2., “Futhermore, the data of the well with precipitation or so on at anytime after the mix of the substance and the cell should be reject for unsuitable suspension” was changed to “

Furthermore, the data of the well with precipitation or so on at anytime after the mix of the

substance and the cell, especially after 72 hr incubation, should be rejected for unsuitable suspension”.

The address of the author, “2-12-1, Fukuura, Kanazawa-ku, Yokohama-shi, 236-8643 Japan” was changed to “2-2-1, Hayabuchi, Tsuzuki-ku, Yokohama-shi, 224-8558, Japan”.

(8)The revision from ver.2.12 from ver.3.1:

The version and the date were renewed.

In 1., “and United States Environmental Protection Agency (EPA)” was deleted.

In 3.1, “(e.g. When the cell culture starts at passage number of 435 and is maintained by two passages per one week, it should be used within passage number of 458)” was added.

In 3.4., “The manufacturer and so on of the reagent is shown in the Table 1.” was changed to “The example of the manufacturer and so on of the reagent is shown in the Table 1.”.

In 3.7.1., “The solubility of the substance in the medium should be confirmed in advance.” was changed to “The solubility of each substance in the medium should be confirmed in advance, using the procedure shown in Fig. 1.”.

In figure 1, “Select the solvent on the basis of the examination of the solubility in DMSO or Ethanol.” was changed to “Examine which solvent is more soluble, DMSO or Ethanol, and select appropriate solvent.”

In 4.6., “For EPA standard, the test substance is judged as negative (Category 4) when the IC50 is higher than that of triethanolamine, and is judged as positive (Category 1-3) when the IC50 is lower than or equal to that of triethanolamine (see table 3 and 4 of appendix 1).” was deleted.

In table 1., “The manufacturer and so on of the reagent” was changed to “The example of the manufacturer and so on of the reagent” . Also, lot numbers were deleted from remarks.

In appendix 1, “The same examination that the in vivo evaluation is performed by EPA classification is shown in table 2. The in vivo results are discriminated between category 4 and others (category 1-3) of the EPA standard. Triethanolamine is classified as category 3. Therefore, when the IC50 of future test substance is higher than that of triethanolamine, it should be evaluated as a non irritant. If the IC50 of the test substance is lower than or equal to that of triethanolamine, it should be evaluated as an irritant.” was deleted.

Table 2 and Table 4 of appendix 1 were deleted and table number was moved up.

(8)The revision from ver.3.1 to ver.3.2:

In 4.7. (4), “Treatment of IC50 expressed with inequality sign is performed in the same manner.” was added.

(9)The revision from ver.3.2 to ver.3.3:

In “1.Purpose”, “as a bottom up approach (Scott et al., 2010)” was added.

In “2. The principle of SIRC cytotoxicity test”, “The relative control was additionally used to obtain the invariable results (Ohno, 2004).” was added.

In triethanolamine of “3.4. Culture medium and reagent”, “It should be used that of Purity $\geq$ 98.0%.” was added.

In 3.7.2., overlap with 3.7.1 was removed. "The test substance is solved or uniformly suspended with medium at a concentration of 10,000 µg/mL (1 w/v%). It is solved or uniformly suspended using vortex mixer, waterbath and sonicator when it finds necessary. DMSO or ethanol is used for solving or suspending if needed. The concentration of DMSO or ethanol is 10,000µg/mL (1 w/v%) in the initial substance solution. The solvent selection is medium, DMSO in medium and ethanol in medium in order. In addition, the concentration of the substance is decreased to 5,000 µg/mL (0.5 w/v%) for suspending when it finds necessary. The substance that is not suspended homogeneously is judged as an inapplicable substance for this test." was deleted. And, "The test substance is solved or uniformly suspended at appropriate concentration by the procedure described in 3.7.1." was added.

In "4.6.Evaluation", "on the basis of the in vivo data by Ohno et al. " was deleted.

In "5.References" and the text, the reference numbers were added. The format of reference was changed. The references of the following paper were added.

Ohno, Y., ATLA, 32, Supplement 1, 643-655, 2004.

Scott, L. et al., Toxicol. *In Vitro*, 24(1), 1-9 (2010).

The reference of the following paper was deleted.

Ohno, Y. et al., *In Vitro Toxicol.*, 7, 89 (1994)

Table 2 of appendix 1 and the related sentences were deleted.

In reference of appendix 1, the format of references was changed.

(10)The revision from ver.3.3 to ver.3.4:

In 2., "The application time of 72hr is set with consideration of time (within 72hr) from ocular instillation to maximal eye irritation for general chemicals except acid or alkali, on the basis of in vivo data at the previous research project<sup>9)</sup>." was added.

In 2., "The SIRC cytotoxicity test procedure is based on the measurement of viable cells stained by crystal violet." was changed to "The SIRC cytotoxicity test is based on the measurement of viable cells stained by crystal violet, which penetrates via cell membrane and stains biological macromolecules."

In 2., "On the other hand, the disadvantage of this method is to be confined to test substances which are solved or uniformly suspended in the medium." was added.

In 4.5., "If the multiple concentrations showing 50% cell viability were obtained from one substance, the lowest IC<sub>50</sub> should be adopted." was added.

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Table 1. The example of the manufacturer and so on of the reagent

No.	Reagent or Medium	Manufacturer	Catalog number	Remarks
1	MEM (Minimum Essential Medium)	GIBCO	Code#: 11095	
2	Fetal Bovine Serum	GIBCO	REF#: 26140-079	
3	Penicillin-Streptomycin-Amphotericin (100x)	GIBCO/BRL	REF#: 15240-062	
4	Phosphate-Buffered Saline (PBS(-))	Nissui	Code#: 05913	
5	0.25% (w/v) Trypsin (1mmol/L EDTA·4Na)	Wako	Cat#: 209-16941	
6	Dimethyl sulfoxide (DMSO)	Kanto	Cat# 2950-1B	
7	Ethanol	Wako	Cat#: 057-00456	
8	Crystal violet	Wako	Cat#: 031-04852	
9	Methanol	Wako	Cat#: 131-01826	
10	Sodium Dodecyl Sulfate	Wako	Cat#: 191-07145	
11	Triethanolamine	Kanto	Cat#: 40268-00	

The products of the same specification are used for the reagents of No.1, 2, 3, 10 and 11. The equivalents regardless of manufacturer are acceptable for other reagents.

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PBS:PBS(-) of 200  $\mu$ L, NC:Medium, 10,000  $\mu$ g/mL DMSO-Medium solution or 10,000  $\mu$ g/mL ethanol-Medium solution of 100  $\mu$ L, S:Dilution series of the test substances by a common ratio of two (100  $\mu$ L), R:Dilution series of the relative control by a common ratio of two (100  $\mu$ L), P: Dilution series of the positive control by a common ratio of two (100  $\mu$ L).


The dilution series of the test substance should be made using medium, 10,000  $\mu$ g/mL DMSO-Medium solution or 10,000  $\mu$ g/mL ethanol-Medium solution. The dilution series of positive control and relative control should be made using medium.

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	1	2	3	4	5	6	7	8	9	10	11	12
A												
B		■										
C												
D												
E												
F												
G												
H												

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Figure 3. Addition of cell suspension

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## Appendix 1 The reason for selecting triethanolamine as a reference control

Triethanolamine was selected as a relative control substance of the SIRC cytotoxicity test for discriminating between non irritant (NI) (= not classified) and irritant in the Globally Harmonized System of Classification and Labeling of Chemicals (GHS). It is one of the substances used in the previous validation study that was performed by a Research Grant for Health Sciences, MHW and Japanese Cosmetic Industry Association, and reported by Ohno et al.<sup>1)</sup> and Tani et al.<sup>2)</sup> Also, it has the ready availability as a commercial reagent, the solvability in the medium, and the useful cytotoxicity and in vivo data for evaluating the eye irritancy.

For selecting the relative control substance, the accuracy and others for discriminating between non irritant and irritant was checked using every available substances in the previous validation with the intention of selecting each substance as a relative control substance. As a result of the comparative study using GHS classification, the triethanolamine was selected as a relative control substance because of the relatively small false negative, the highest accuracy except for substances which could not be obtained the cytotoxicity (10000< etc.), and the clarity about the category (for example, alcohol) of false negative substances, as shown in table 1. The identification of the group of the false negative substances is important to avoid the risk of eye irritation.

### References

- 1) Ohno, Y. et al., *Toxicology in Vitro* 13, 73-98 (1999)
- 2) Tani, N. et al., *Toxicology in Vitro* 13, 175-187(1999)

Table 1 The correlative evaluation on the basis of the previous validation study data for selecting relative control substance – GHS classification –

Substances	In vivo Evaluation by GHS <sup>§</sup>	The SIRC cytotoxicity : IC <sub>50</sub> (ug/mL) <sup>#</sup>	The values for ranking	True Negative	False Negative	True Positive	False Positive	Accuracy (%)
Polyethylene glycol 400	NI	35300<	35300	1	0	27	6	82
Silicic anhydride	NI	14800<	14800	2	0	26	6	82
Glycerin	NI	11600	11600	3	0	25	6	82
Isotonic sodium chloride solution	NI	10000<	10000	4	1	24	5	82
Ethanol	1or2A	10000<	10000	4	1	24	5	82
Isopropyl myristate	NI	9330<	9330	5	1	24	4	85
Butanol*	1or2A	8880<	8880	5	2	23	4	82
Triethanolamine	NI	2090	2090	6	2	23	3	85
Lactic acid	1	1230	1230	6	3	22	3	82
Benzyl alcohol	1or2A	1190	1190	6	4	21	3	79
Polyoxyethylene sorbitan monooleate (20E.O.)	NI	963	963	7	4	21	2	82
Sodium salicylate	1or2A	952	952	7	5	20	2	79
Glycolic acid*	1or2A	868	868	7	6	19	2	76
Acetic acid*	1or2A	721	721	7	7	18	2	74
Diisopropanolamine*	1, 2Aor2B	699	699	7	8	17	2	71
2-Ethylhexyl p-dimethylamino benzoate	NI	474	474	8	8	17	1	74
Calcium thioglycolate	1	392	392	8	9	16	1	71
Acid red 92	1or2A	297	297	8	10	15	1	68
Sucrose fatty acid ester	1or2A	286	286	8	11	14	1	65
m-Phenylenediamine	1or2A	218	218	8	12	13	1	62
Methyl p-hydroxybenzoate	NI	207	207	9	12	13	0	65
Di (2-ethylhexyl) sodium sulfosuccinate*	1or2A	181	181	9	13	12	0	62
Sodium lauryl sulfate*	1or2A	168	168	9	14	11	0	59
Sodium hydrogenated tallow L-glutamate*	1or2A	140	140	9	15	10	0	56
Potassium laurate*	1or2A	120 (Data from 4 labs)	120	9	16	9	0	53
Chlorhexidine gluconate (20% solution)*	1or2A	67.6	67.6	9	17	8	0	50
Polyoxyethylene octylphenylether (10 E.O.)*	1or2A	38.4	38.4	9	18	7	0	47
Distearyldimethylammonium chloride	1	37.8	37.8	9	19	6	0	44
Benzalkonium chloride*	1or2A	19.0	19	9	20	5	0	41
Domiphen bromide*	1or2A	12.1	12.1	9	21	4	0	38
Monoethanolamine*	1or2A	9.62	9.62	9	22	3	0	35
Cetyltrimethylammonium bromide*	1or2A	2.59 (Data from 4 labs)	2.59	9	23	2	0	32
Cetylpyridinium chloride*	1	1.67	1.67	9	24	1	0	29
Stearyltrimethylammonium chloride*	1	1.58	1.58	9	25	0	0	26

§ : The Draize eye test results couldn't always discriminate between 1 and 2 of GHS classification for no observation data on day 21. The observation was performed to day 14.

# : The data of SIRC cytotoxicity test were the mean of IC<sub>50</sub> (ug/mL) from more than 5 laboratory except for one part.

\* : The in vivo results of as is application was predicted from the data of 10% concentration.

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Appendix 2 The basis for the set IC<sub>50</sub> range of triethanolamine as a reference control

The IC<sub>50</sub> range of triethanolamine, 1,000-2,500 µg/mL is based on the mean  $\pm$  2 x standard deviation from the data (n=144).