

新規試験法提案書

眼刺激性試験代替法
SkinEthic™ HCE TTT法

令和6年2月

国立医薬品食品衛生研究所

新規試験法提案書

令和 6 年 2 月 26 日

No. 2023-02

眼刺激性試験代替法SkinEthic™ HCE TTT法に関する提案

令和 5 年 12 月 12 日に国立医薬品食品衛生研究所にて開催された新規試験法評価会議（通称：JaCVAM 評価会議）において以下の提案がなされた。

提案内容： 本試験法は、技術移転性、施設内再現性、施設間再現性の高い試験法である。予測性に関しては、適用範囲に留意すれば、UN GHS 区分 1、区分 2 への分類および区分に該当しない場合の判定を可能とすると考えられた。ただし、UN GHS 区分 2A と 2B を区別することはできない。また、固体の UN GHS 区分 2 を区分に該当しないと判断する偽陰性率が 28.9%と高く、30%以下という OECD が定めた基準に近い値であることから、固体の UN GHS 区分 2 の分類は慎重に評価されるべきである。

この提案書は、眼刺激性試験資料編纂委員会によりまとめられた文書を用いて、JaCVAM 評価会議が評価および検討した結果、その有用性が確認されたことから作成された。

以上の理由により、行政当局の安全性評価方法として眼刺激性試験代替法 SkinEthic™ HCE TTT 法の使用を提案するものである。



西川秋佳

JaCVAM 評価会議 議長



平林容子

JaCVAM 運営委員会 委員長

JaCVAM 評価会議

- 西川 秋佳 (国立医薬品食品衛生研究所 安全性生物試験研究センター 病理部／
名古屋徳洲会総合病院) : 座長
- 小島 幸一 (一般財団法人 食品薬品安全センター)
- 中村 りこ (独立行政法人 製品評価技術基盤機構)
- 西村 次平 (独立行政法人 医薬品医療機器総合機構)
- 平林 容子 (国立医薬品食品衛生研究所 安全性生物試験研究センター)
- 松本 一彦 (名古屋市立大学大学院)

任期：令和4年4月1日～令和6年3月31日

JaCVAM 運営委員会

- 平林容子 (国立医薬品食品衛生研究所 安全性生物試験研究センター) : 委員長
石井孝司 (国立感染症研究所)
稲角嘉彦 (厚生労働省 医薬・生活衛生局 医薬品審査管理課 化学物質安全対策室)
小川久美子 (国立医薬品食品衛生研究所 安全性生物試験研究センター 病理部)
諫田泰成 (国立医薬品食品衛生研究所 安全性生物試験研究センター 薬理部)
北嶋 聡 (国立医薬品食品衛生研究所 安全性生物試験研究センター 毒性部)
杉山圭一 (国立医薬品食品衛生研究所 安全性生物試験研究センター 変異遺伝部)
高橋祐次 (国立医薬品食品衛生研究所 安全性生物試験研究センター 毒性部 動物管理室)
束野正明 (厚生労働省 医薬・生活衛生局 医薬品審査管理課 化学物質安全対策室)
林 亜紀子 (厚生労働省 医薬・生活衛生局 医薬品審査管理課)
本間正充 (国立医薬品食品衛生研究所)
真木一茂 (独立行政法人 医薬品医療機器総合機構)
増村健一 (国立医薬品食品衛生研究所 安全性生物試験研究センター 安全性予測評価部)
横田雅彦 (独立行政法人 医薬品医療機器総合機構)
足利太可雄 (国立医薬品食品衛生研究所 安全性生物試験研究センター 安全性予測評価部
第二室) : 事務局

JaCVAM statement on the SkinEthic™ HCE TTT, an alternative method for evaluating eye irritation

At a meeting held on 12 December, 2023 at National Institute of Health Sciences (NIHS) in Tokyo, Japan, the Japanese Center for the Validation of Alternative Methods (JaCVAM) Regulatory Acceptance Board unanimously endorsed the following statement:


Proposal: This test method has high technological transferability and within- and between-laboratory reproducibility. Regarding predictive capacity, we believed that consideration of the applicability domain would enable the classification of a substance as United Nations Globally Harmonized System of Classification and Labeling of Chemicals (UN GHS) category 1 or 2 or determine no category. However, it was impossible to distinguish between UN GHS categories 2A and 2B. Moreover, regarding a solid, as the false negative rate for discriminating UN GHS category 2 and no category is high (28.9%; close to the standard set by the OECD of $\leq 30\%$), the classification of UN GHS category 2 for solids should be performed with caution.

This statement was released following a review prepared by the eye irritation test JaCVAM Editorial Committee to acknowledge that the results of the review and study by the JaCVAM Regulatory Acceptance Board have confirmed the usefulness of this assay.

Based on the above, we proposed the SkinEthic™ HCE TTT method as a useful means for assessing eye irritation potential during safety assessments by regulatory agencies.



Nishikawa Akiyoshi
Chairperson,
JaCVAM Regulatory Acceptance Board.



Hirabayashi Yoko
Chairperson,
JaCVAM Steering Committee.

February 26, 2024

The JaCVAM Regulatory Acceptance Board was established by the JaCVAM Steering Committee, and is composed of nominees from the industry and academia.

This statement was endorsed by the following members of the JaCVAM Regulatory Acceptance Board:

Nishikawa Akiyoshi (Division of Pathology, Center for Biological Safety and Research:CBRS,
NIHS / Nagoya Tokushukai General Hospital) : Chairperson

Hirabayashi Yoko (CBRS, NIHS)

Kojima Koichi (Food and Drug Safety Center)

Matsumoto Kazuhiko (Nagoya City University)

Nakamura Ruriko (National Institute of Technology and Evaluation)

Nishimura Jihei (Pharmaceuticals and Medical Devices Agency)

Term: From 1st April 2022 to 31st March 2024

This statement was endorsed by the following members of the JaCVAM steering Committee after receiving the report from JaCVAM Regulatory Acceptance Board:

Hirabayashi Yoko (CBSR, NIHS): Chairperson

Hayashi Akiko (Ministry of Health, Labour and Welfare)

Honma Masamitsu (NIHS)

Inazumi Yoshihiko (Ministry of Health, Labour and Welfare)

Ishii Koji (National Institute of Infectious Diseases)

Kanda Yasunari (Division of Pharmacology, CBSR, NIHS)

Kitajima Satoshi (Division of Cellular and Molecular Toxicology, CBSR, NIHS)

Maki Kazushige (Pharmaceuticals and Medical Devices Agency)

Masumura Kenichi (Division of Risk Assessment, CBSR, NIHS)

Ogawa Kumiko (Division of Pathology, CBSR, NIHS)

Sugiyama Keiichi (Division of Genetics and Mutagenesis, CBSR, NIHS)

Taquahashi Yuhji (Animal Management Section of the Division of Toxicology, CBSR, NIHS)

Tsukano Masaaki (Ministry of Health, Labour and Welfare)

Yokota Masahiko (Pharmaceuticals and Medical Devices Agency)

Ashikaga Takao (Division of Risk Assessment, CBSR, NIHS): Secretary

眼刺激性試験代替法SkinEthic™ HCE TTT法に関する提案

添付資料一覧

添付資料 1

評価会議報告書 ----- 1

添付資料 2

評価報告書 ----- 5

添付資料 3

OECD GUIDELINE FOR TESTING OF CHEMICALS 492B
RECONSTRUCTED HUMAN CORNEA-LIKE EPITHELIUM (RHCE) TEST
METHOD FOR EYE HAZARD IDENTIFICATION-----21

添付資料 1

評価会議報告書

眼刺激性試験代替法 SkinEthic™ HCE TTT 法

JaCVAM 評価会議

令和 5 年(2023 年)12 月 12 日

JaCVAM 評価会議

- 西川 秋佳 (国立医薬品食品衛生研究所 安全性生物試験研究センター 病理部／
名古屋徳洲会総合病院) : 座長
- 小島 幸一 (一般財団法人 食品薬品安全センター)
- 中村 りこ (独立行政法人 製品評価技術基盤機構)
- 西村 次平 (独立行政法人 医薬品医療機器総合機構)
- 平林 容子 (国立医薬品食品衛生研究所 安全性生物試験研究センター)
- 松本 一彦 (名古屋市立大学大学院)

任期：令和4年4月1日～令和6年3月31日

JaCVAM 評価会議は、眼刺激性試験資料編纂委員会により作成された「SkinEthic™ HCE TTT 法評価報告書」¹⁾をもとに本試験法の科学的妥当性、社会的および行政的な受け入れ性について検討した。

1. 試験法の定義および科学的妥当性

名称：SkinEthic™ Human Corneal Epithelium (HCE) Time-to-Toxicity (TTT)法

代替する対象毒性試験：Draize 眼刺激性試験²⁾

科学的妥当性：

当該試験法は、再構築ヒト角膜様上皮モデル(Reconstructed human Cornea-like Epithelium: RhCE)である SkinEthic™ HCE 用い、細胞毒性を指標として眼刺激性を評価する試験法である。経済協力開発機構(Organisation for Economic Co-operation and Development: OECD)試験法ガイドライン(Test Guideline: TG)492B に記載されている³⁾。

液体または固体の化学物質毎に異なる濃度および培養時間にて処理することにより生じる細胞生存率を用い、国際連合化学品の分類および表示に関する世界調和システム(United Nations Globally Harmonized System of Classification and Labelling of Chemicals: UN GHS)による UN GHS 区分 1、区分 2 への分類および区分に該当しない場合の判定ができる方法として開発されており、科学的には妥当である。

2. 目的とする物質又は製品の毒性を評価する試験法としての、社会的受け入れ性および行政上の利用の可能性

社会的受け入れ性：

本試験法は RhCE に対する化学物質の細胞毒性を指標に用いて眼刺激性を評価する試験法であり、生きた動物を用いないという点で、3Rs の精神に合致している。また、SkinEthic™ HCE TTT 法に用いるキットの入手は容易であり、短時間で実施でき、特殊な機材や試薬を必要とせず、必要な手技も複雑なものでない。したがって、入手したキットの品質基準が許容範囲にあり、かつ実施する試験施設の技術習得がガイドラインの熟達度確認物質で確認できており、基本的な細胞培養の技術と設備を有する施設であれば実施可能である。技術移転性は高い上に、原料だけでなく、製剤においても UN GHS 区分 1、区分 2 への分類および区分に該当しない場合の判定に用いることができる。以上より、本試験法の社会的受け入れ性は高い。

行政上の利用性：

本試験法は、技術移転性、施設内再現性、施設間再現性の高い試験法である。予測性に関しては、適用範囲に留意すれば、UN GHS 区分 1、区分 2 への分類および区分に該当しない場合の判定を可能とすると考えられた。ただし、UN GHS 区分 2A と 2B を区別することはできない。また、固体の UN GHS 区分 2 を区分に該当しないと判断する偽陰性率が 28.9% と高く、30%以下という OECD が定めた基準に近い値であることから、固体の UN GHS 区分 2 の分類は慎重に評価されるべきである。

参考文献（最終確認日：2023年9月20日）

- 1) JaCVAM 眼刺激性試験資料編纂委員会：眼刺激性試験代替法SkinEthic™HCE TTT法評価報告書（2023年3月20日）
- 2) OECD (2023) OECD Guidelines for the Testing of Chemicals No. 405. Acute Eye Irritation/Corrosion, Organisation for Economic Cooperation and Development, Paris.
Available at: <https://www.oecd-ilibrary.org/docserver/9789264185333-en.pdf?expires=1695282151&id=id&accname=guest&checksum=3130BC4DE84C35C7D73325F595BDA84B>
- 3) OECD (2022) OECD Guidelines for the testing of Chemicals No. 492B. Reconstructed human Cornea-like Epithelium (RhCE) test method for Eye Hazard Identification, Organisation for Economic Cooperation and Development, Paris.
Available at: https://www.oecd-ilibrary.org/environment/test-no-492b-reconstructed-human-cornea-like-epithelium-rhce-test-method-for-eye-hazard-identification_0d603916-en

添付資料 2

評価報告書

眼刺激性試験代替法 SkinEthic™ HCE TTT 法

眼刺激性試験資料編纂委員会

令和5年（2023年）12月12日

眼刺激性試験資料編纂委員会

山本直樹（委員長：藤田医科大学）

小島 肇（国立医薬品食品衛生研究所）

佐々木正治（アレクシオンファーマ合同会社）

竹内小苗（P&G イノベーション合同会社）

波多野浩太（ホーユー株式会社）

原 範子（藤田医科大学）

山下晴洋（アステラス製薬株式会社）

略語

CAS:	Chemical Abstracts Services
GHS:	Globally Harmonized System of Classification and Labelling of Chemicals
HCE:	Human Corneal Epithelium
HPLC:	High Performance Liquid Chromatography
IATA :	Integrated Approaches to Testing and Assessment
JaCVAM:	Japanese Center for the Validation of Alternative Methods
MTT:	3- (4,5-Dimethyl-2-thiazolyl)-2,5-diphenyltetrazolium Bromide
OECD:	Organization for Economic Co-operation and Development
PBS:	Phosphate-Buffered Saline
PRP:	Peer Review Panel
RhCE:	Reconstructed human Cornea-like Epithelium
TG:	Test Guideline
TTT:	Time-to-Toxicity
UN:	United Nations
UPLC:	Ultra Performance Liquid Chromatography
VRM:	Validated Reference Method

要旨

SkinEthic™ Human Corneal Epithelium (HCE) Time-to-Toxicity (TTT)法は、L'Oréal 社によって開発され、再構築ヒト角膜様上皮モデル(Reconstructed human Cornea-like Epithelium: RhCE)である SkinEthic™ HCE 用いて、液体および固体の化学物質における国際連合化学品の分類および表示に関する世界調和システム(United Nations Globally Harmonized System of Classification and Labelling of Chemicals: UN GHS)による UN GHS 区分 1、UN GHS 区分 2 への分類および区分に該当しない場合の判定ができる方法として開発された眼刺激性評価試験法である。

SkinEthic™ HCE TTT 法は、バリデーション研究の後、科学的妥当性を評価するために第三者評価委員会(Peer Review Panel : PRP)による検証がなされた。その後、経済協力開発機構(Organization for Economic Co-operation and Development: OECD)による専門家会議での議論を経て、UN GHS 区分 1、UN GHS 区分 2 への分類および区分に該当しない場合の判定が可能な試験法と判断され、OECD 試験法ガイドライン(Test Guideline : TG)492B として 2022 年 6 月に採択された。

本バリデーション研究は、本試験法に用いられる RhCE である SkinEthic™ HCE を販売している企業の親会社である L'Oréal 社がスポンサーとなり行われた。そのため、バリデーション実行委員会の独立性について、PRP で慎重に議論された。その結果、PRP はバリデーション報告書の記載に不備は多いが、バリデーションの結果は受け入れると結論している。いずれの施設の記録も十分であり、データも適切であることを PRP が確認していることから、バリデーション結果に大きな問題はないと本委員会も考える。

本試験法の技術移転性は高いと判断された。再現性に関しては、液体被験物質(Time-To-Toxicity Test Method On Liquids : TTL)について、3 施設でのコード化された 20 物質によるバリデーション研究の結果から、施設内再現性は 85-95%であり、施設間再現性は 90%であった。固体被験物質(Time-To-Toxicity Test Method On Solids: TTS)について、3 施設でのコード化された 20 物質によるバリデーション研究の結果から、施設内再現性は 100%であり、施設間再現性も 100%であったことから、SkinEthic™ HCE TTT 法の再現性は高いと考えられた。また、本試験法の正確度については、UN GHS 区分 1 物質、区分 2 物質、および区分に該当しない物質がそれぞれ 79.2%、69.2%および 74.9%の割合で正しく判定され、本試験法の予測性は OECD 専門家会議が提案した採用基準を満たしていた。以上より、SkinEthic™ HCE TTT 法は、技術移転性、施設内再現性、施設間再現性は高い試験法である。ただし、UN GHS 区分 2 の固体物質を区分に該当しない物質と判定する偽陰性率が 28.9%と高く、30%以下という基準に近い値であることから、安全性評価上懸念を持つとの見解で一致した。

以上より、SkinEthic™ HCE TTT 法は、OECD TG492B に準拠して実施した場合、UN GHS 区分 1、UN GHS 区分 2 への分類および区分に該当しない場合の判定を可能とする試験法である。ただし、固体物質の UN GHS 区分 2 の分類は慎重に評価されるべきであると本資料編纂委員会は考える。

1. 緒言

OECD は TG として、既に複数の眼刺激性試験代替法を採択している¹⁾。UN GHS 区分 1 の分類および UN GHS 区分に該当しない場合を判定試験法としてウシ摘出角膜の混濁および透過性試験法 (Bovine Corneal Opacity and Permeability test method: BCOP 法、TG437)、ニワトリ摘出眼球を用いた眼刺激性試験 (Isolated Chicken Eye Test: ICE 法、TG438)、*in vitro* 短時間暴露法 (Short Time Exposure: STE 法、TG491)、*in vitro* マクロモレキュラ試験法 (TG496) があり、UN GHS 区分 1 を検出する試験法としてフルオレセイン漏出試験法 (TG460)、UN GHS 区分に該当しない場合を判定する試験法として再構築ヒト角膜様上皮モデル法 (Reconstructed human Cornea-like Epithelium: RhCE 法、TG492)²⁾ および Vitrigel-Eye Irritancy Test (Vitrigel-EIT, TG494) がある。ただし、いずれの試験法においても、UN GHS 区分 2 の化学物質の検出をすることはできず、*in vivo* ドレイズ試験の代替法ではない。

この度、UN GHS 区分 1、UN GHS 区分 2 への分類および区分に該当しない場合の判定が可能となる試験法として L'Oréal 社によって開発された SkinEthic™ Human Corneal Epithelium (HCE) Time-to-Toxicity (TTT) 法は、RhCE である SkinEthic™ HCE を用いて液体および固体に対する化学物質の眼刺激性の評価を行う試験法である^{3,4)}。SkinEthic™ HCE TTT 法は、L'Oréal 社がスポンサーとなり、L'Oréal R&I と Adriaens Consulting BVBA の主導で実施されたバリデーション研究の後、科学的妥当性を評価するために第三者評価委員会 (Peer Review Panel: PRP) による検証がなされた。その後、OECD による専門家会議での議論を経て、UN GHS 区分 1、UN GHS 区分 2 への分類および区分に該当しない場合の判定が可能となる試験法と判断され、OECD TG492B として 2022 年 6 月に採択された⁵⁾。

本報告書は、バリデーション研究論文⁶⁾、第三者評価報告書⁷⁾、その他関連論文などをもとに本試験法の概要を説明し、JaCVAM 眼刺激性試験資料編纂委員会の意見をまとめたものである。

2. 試験法の位置づけ

SkinEthic™ HCE TTT 法は、化学物質 (単一物質および混合物) の眼刺激性を UN GHS において UN GHS 区分 1、UN GHS 区分 2 への分類および区分に該当しない場合の判定に用いることができる試験法である。

3. 原理

本試験法は、SkinEthic™ HCE を用いて細胞毒性を指標として眼刺激性を評価する試験法である。角膜混濁などの細胞傷害は種々の作用機序に伴って発生する可能性があるが、細胞毒性は *in vivo* で見られる化学物質の眼刺激性を決定するうえで重要な役割を果たす。また、化学物質の眼刺激性は主に初期損傷の程度によって決定され、細胞死の程度とその後の反応および最終的な結果の程度と相関する。

4. 試験手順

SkinEthic™ HCE TTT 法の手順を以下に示す^{8,9)}。SkinEthic™ HCE は EPISKIN 社(フランス)および我が国では株式会社ニコダームリサーチより購入できる。

4-1. RhCE 組織の機能的条件

RhCE 組織の機能的条件は、以下のとおりである。

陰性対照吸光度：陰性対照において $1.0 < OD \leq 2.5$

(使用時の基準であり、製造者の出荷基準は異なる場合もある)

バリア機能：50 μ L ドデシル硫酸ナトリウム(Sodium Lauryl Sulphate : SDS)の30分間曝露において $1.0 \leq IC_{50}(\text{mg/mL}) \leq 3.2$

(IC_{50} ：細胞生存率を50%低下させるのに必要な濃度)

形態：少なくとも4層の上皮細胞があり、角化していない

再現性：陽性および陰性対照の結果が背景データをもとに設定した許容範囲内にある

これらは SkinEthic™ HCE 製造者の出荷基準として採用される。一方、SkinEthic™ HCE の使用者は毎回、生存率と再現性を確認する必要がある。

4-2. 被験物質の適用

37°C 以下においてピペットで扱えるものは液体として、それ以外は固体として試験を行う。RhCE 組織は、一晚、標準培養条件(37 \pm 2°C、5 \pm 1% CO₂、 \geq 95% 湿度)で培養し、更に新たな培地を用いて標準培養条件で30分間以上培養する。

液体物質の場合(SkinEthic™ HCE Time-to-Toxicity for Liquids chemicals: TTL)、評価対象物質そのもの、および蒸留水で20%(w/v)に希釈したものをを用いる。SkinEthic™ HCE に Phosphate-Buffered Saline (PBS) 10 \pm 2 μ L を適用し被験物質 80 \pm 2 μ L を加える。評価対象物質そのものは室温で5 \pm 0.25分間、20%(w/v)希釈調整液は室温で16 \pm 1分間および標準培養条件で120 \pm 2分間培養する。その後、PBS 25mL で洗浄して被験物質を除去し、10 \pm 1分間室温で培地に浸漬する。同時に行う陽性対照には、酢酸メチル(Chemical Abstracts Services : CAS No. 79-20-9)が推奨され、適用・培養・後処理は対象評価物質の条件に準じる。同時に行う陰性対照には PBS が推奨され、80 \pm 2 μ L のみを適用する。培養・後処理は対象評価物質の条件に準じる。

固体物質の場合(SkinEthic™ HCE Time-to-Toxicity for Solids chemicals: TTS)、評価対象物質そのものを用いる。被験物質 80 \pm 2 mg を均一に適用した後に蒸留水 80 \pm 2 μ L を加え、標準培養条件で30 \pm 2分間および120 \pm 5分間培養する。その後、PBS 25 mL で洗浄して被験物質を除去し、30 \pm 2分間室温で培地に浸漬する。同時に行う陽性対照として、乳酸1%(w/v)水溶液(CAS No. 50-21-5)が推奨され、80 \pm 2 μ L を適用した後に蒸留水 80 \pm 2 μ L を加える。同時に行う陰性対照には PBS が推奨され、80 \pm 2 μ L のみを適用する。陽性・陰性対照とも培養・後処理は対象評価物質の条件に準じる。

それぞれの曝露条件で少なくとも2つのRhCE組織を用いる。

4-3. 細胞生存率の算出

細胞生存率算出には 3-(4,5-di-methylthiazol-2-yl)-2,5-diphenyltetrazolium bromide (MTT) 還元法を用いる。培地を除去し、1 mg/mL MTT 溶液 0.3 mL 中で 180±15 分間標準培養条件で反応させ、その後 1.5 mL イソプロパノール(または同様な溶媒)で青色のホルマザンを抽出する。液体被験物質の場合は、RhCE 組織の上部・底部両方から抽出する。固体被験物質および無色透明でない液体被験物質の場合は、組織に残存する被験物質の混入を最小限に抑えるため、ホルマザンの抽出は RhCE 組織の底部のみから行う。液体被験物質でも、洗浄が困難な場合には、底部のみから抽出を行う。同時対照物質に対しての抽出方法は被験物質と同様に行う。抽出したホルマザンの定量は OD_{570nm} 測定または High Performance Liquid Chromatography (HPLC)/Ultra Performance Liquid Chromatography (UPLC)で行う。

被験物質が MTT 還元物質の場合、あるいはホルマザンと同じような波長 (570nm±30nm) に吸収を持つ着色物質の場合、細胞生存率の補正を行う必要がある。その手順については TG 本文の説明 (33-40 項) および Annex 3 のフローチャートを参照する⁵⁾。

4-4. 試験成立の承認基準

以下の条件をすべて満たした場合、試験の成立を承認する。

- 1) 陰性対照の平均 OD が $1.0 < OD \leq 2.5$ であること。
- 2) 液体被験物質では、陽性対照物質の平均細胞生存率が 5 分間の曝露時間で 50%以下、16 分および 120 分間の曝露時間で 50%を超えること、固体被験物質では、陽性対照物質は平均細胞生存率が 30 分間の曝露時間で 40%を超えること、120 分間の曝露時間で $20\% < \% \text{生存率} \leq 70\%$ であること。いずれも背景データの範囲内にあること。
- 3) 被験物質および陰性・陽性対照のそれぞれにおいて 2 つの RhCE 組織の細胞生存率の差が 20%未満であること。3 つの RhCE 組織用いた場合は細胞生存率の標準偏差が 18%以下であること。

なお、複数の被験物質を同時に試験した場合は、この条件を満たさない被験物質のみ不成立とする。

4-5. 刺激性の判定

液体被験物質では、すべての曝露時間 (5 分、16 分および 120 分) で平均細胞生存率が 50%を超えた場合、被験物質は UN GHS 区分に該当しない、すべての曝露時間で平均細胞生存率が 50%以下の場合 UN GHS 区分 1 と判断される。それら以外の場合は UN GHS 区分 2 と判断される。たとえば、5 分間(原体)と 16 分間(20%希釈液)の曝露条件では平均細胞生存率が 50%を超えていても、120 分間(20%希釈液)で 50%以下である物質は UN GHS 区分 2 となる。

固体被験物質では、平均細胞生存率が 30 分間(原体)の曝露時間で 40%を超え、かつ 120 分間(原体)の曝露時間で 60%を超えた場合、被験物質は UN GHS 区分に該当しないと判断され、平均細胞生存率が 30 分間(原体)の曝露時間で 40%以下、かつ 120 分間(原体)の曝露時間で 60%以下の場合は UN GHS 区分 1 と判断される。それら以外の場合は UN GHS 区分 2 と判断される。

なお、平均細胞生存率がカットオフ値近辺の場合(±5%)は 2 回目の試験を検討する。1 回目と 2 回目の試験で結果が一致しない場合は 3 回目の試験を検討する。

5. バリデーション研究

5-1. バリデーション研究の透明性

本バリデーション研究は、本試験法に用いられる RhCE である SkinEthic™を販売している企業の親会社である L'Oréal 社がスポンサーとなり行われた。そのため、バリデーション実行委員会の独立性が PRP でも慎重に議論された⁷⁾。バリデーション実行委員会の中の第三者機関としては、コンサルティング企業である Adriaens Consulting BVBA のみが発関し、コード化やデータ解析を実施しているが、被験物質の選択は L'Oréal 社が実施している。さらに、統計学者の報告書もバリデーションの採用基準も計画にはなく、バリデーション報告書にも不明確な記載が多かった。これらの事実から、PRP は OECD GD 34 で推奨されている試験マネジメントや被験物質の選択に関する独立性については不明確であるとしている。

このような状況を PRP が慎重に検討して明らかにし、以下の 15 項目について確認した⁷⁾。その結果、PRP はバリデーション報告書の記載に不備は多いが、バリデーションの結果は受け入れると結論している。

- 1) バリデーション研究と試験法の目的(Study objective and test method purpose)
- 2) 既存の試験法との比較による必要性和有用性(Need and benefits in comparison to existing test methods)
- 3) 性能標準(OECD GD216)¹⁰⁾に記載された要件に叶う試験法との比較(Comparison of the test method with the essential test method components as described in the Performance Standards)
- 4) 生物学のおよび作用機序的な関連性 (Biological and mechanistic relevance)
- 5) 試験法プロトコル (Test method protocol)
- 6) バリデーション研究の管理と実行の適正 (Appropriateness of the validation study management and conduct)
- 7) バリデーション研究の目的に沿った適切な被験物質の選択 (Adequacy of chemicals selection for the validation study objective)
- 8) 関連性評価に使われた参照データの質 (Quality of the reference data used for the evaluation of relevance)
- 9) 経験のない施設へのトレーニングと技術移転 (Training of naive laboratories and transferability)

- 10) データの信頼性保証 (Use of quality assurance system (s) during data generation)
- 11) 施設内および施設間再現性 (Within- and between-laboratory reproducibility)
- 12) 予測性 (Predictive capacity)
- 13) 適用範囲と限界 (Applicability domain and limitations)
- 14) データと文書の完成度 (Completeness of data and documentation)
- 15) 他の関連事項 (Other relevant aspects)

なお、予測性に関しては OECD の専門家会議で定められた採用基準をもとに評価された。

本委員会では、バリデーション研究を遂行する場合には、開発者でなく、第三者機関が主導し、透明性や利益相反に考慮するべきであると考えます。ただし、バリデーション報告書の記載に不備は多いが、いずれの施設の記録も十分であり、データも適切であることを PRP が確認しており、バリデーション結果に大きな問題はないと本委員会も考えた。

5-2. トレーニングと技術移転

SkinEthic™ HCE TTT 法のバリデーション研究にはリードラボの EPISKIN 社と本試験法の経験がない 2 施設 (EUROSAFE 社、VITROSCREEN 社) を加えた合計 3 施設が参加した。すべての施設が開発者から 3~4 日間のトレーニングを受けた。

技術移転性として、コード化された 3 物質、陰性対照物質および陽性対照物質を用い、少なくとも 2 回の実験が実施された。液体被験物質の場合 (TTL) において、1 回目で 1 物質の 5 分間処理の判定結果が 2 施設で食い違ったが、細胞毒性 50% の基準前後のバラツキによるものと判断された。しかし、すべての参加施設で 2 回目は正確に予測できていた。固体被験物質の場合 (TTS) では、3 施設の判定結果はすべて同じであった。

以上の結果から、本試験法の技術移転性は高いと判断されている。なお、本試験法を実施する試験施設の技術習得を確かめるための習熟度確認物質の一覧は Appendix 1 (液体は Appendix 1-1、固体は Appendix 1-2) にある。

5-3. 施設内および施設間再現性^{3,4)}

TTL に関しては、コード化された 20 物質 (新規 14 物質を含む) を用いたバリデーション研究が実施され、合計 60 回 (20 物質×3 回) の試験が行われた。施設内再現性は 85、90 および 95% であり、20 物質の施設間再現性は 90% であった。TTS に関しても、コード化された 20 物質 (新規 11 物質を含む) を用いたバリデーション研究において、施設内再現性は 100% であり、施設間再現性も 100% であった。これらの結果から、SkinEthic™ HCE TTT 法の再現性は高いと考えられた。

5-4 試験法の予測性^{1,2,6)}

表 1 に SkinEthic™ HCE TTT 法の正確性を示す。SkinEthic™ HCE TTT 法の正確性は、加重算出した。たとえば、1 種類の被験物質について 3 回の試験結果がある場合で、2 回は Cat.1、1 回は Cat.2 と判定された場合、Cat.1 および Cat.2 の予測はそれぞれ 0.66 (2/3) および 0.33 (1/3) とした。

TTL では、UN GHS 区分 1 物質、UN GHS 区分 2 物質および区分に該当しない物質が、それぞれ、85.4%、79.8% および 79.2% の割合で、また、TTS では、それぞれ、74.7%、55.3% および 71.7% の割合で正しく判定された。TTL と TTS を合わせた TTT の正確度については、UN GHS 区分 1 物質、UN GHS 区分 2 物質および区分に該当しない物質が、それぞれ 79.2%、69.2% および 74.9% の割合で正しく判定され、本試験法の予測性は OECD 専門家会議が承認した採用基準を満たしていた。詳細な結果は Annex 3 および 4 に示している。この基準は SkinEthic™ HCE TTT の開発者の提案を受け、peer review report, 3.12 予測性 (Evaluation criterion 12)¹¹⁾での検討をもとに OECD 専門家会議が承認したものであるが、本委員会でも再度慎重に議論した。特に、UN GHS 区分 2 の固体物質を区分に該当しない物質と判定する偽陰性率が 28.9% と高く、OECD が定めた 30% 以下という採用基準に近い値であることから、安全性評価上懸念を持つとの見解で一致した。

SkinEthic™ HCE TTT 法は、OECD TG492B に準拠して実施した場合、UN GHS 区分 1、UN GHS 区分 2 への分類および区分に該当しない場合の判定を可能とする試験法である。ただし、固体物質の UN GHS 区分 2 の分類は慎重に評価されるべきであるとの意見でまとまった。

表 1 SkinEthic™ HCE TTT 法の正確性*

UN GHS Cat. (# of Chemical)	SkinEthic™ HCE TTL (%)		
	Cat. 1	Cat.2	No Cat.
Cat 1 (21)	85.4	14.6	0.0
Cat. 2 (25)	20.2	79.8	0.0
No Cat. (24)	0.0	20.8	79.2
UN GHS Cat. (# of Chemical)	SkinEthic™ HCE TTS (%)		
	Cat. 1	Cat.2	No Cat.
Cat 1 (29)	74.7	25.3	0.0
Cat. 2 (19)	15.8	55.3	28.9
No Cat. (33)	3.0	25.3	71.7
UN GHS Cat. (# of Chemical)	SkinEthic™ HCE TTT (%)		
	Cat. 1	Cat.2	No Cat.
Cat 1 (50)	79.2	20.8	0
Cat. 2 (44)	18.3	69.2	12.5
No Cat. (57)	1.8	23.4	74.9

*: 値は加重算出した。網掛け箇所が正しく予測できた割合。

of Chemical : 物質数、Cat.1 : 区分 1、Cat.2 : 区分 2、No Cat. : 区分に該当しない

表 2 に SkinEthic™ HCE TTT 法とその他の RhCE 法の UN GHS 区分に該当しない物質に対する正確性を示す。EpiOcular™ Eye Irritation Test (EpiOcular), SkinEthic™ Human Corneal Epithelium Eye Irritation Test (SkinEthic), LabCyte CORNEA-MODEL24 Eye Irritation Test (LabCyte)の正確度、感度および特異度は 80%~88%、95%~100%および 63%~73%であった。SkinEthic™ HCE TTT 法の US GHS 区分に該当しない物質に対する正確度、感度および特異度は、それぞれ、87%、94%および 75%であった。SkinEthic™ HCE TTT 法の正確度、感度、特異度はその他の RhCE 法と同等であった。

表2 UN GHS 区分に該当しない物質に対しての正確性の比較

	正確度 % (# of Chemical)	感度 % (# of Chemical)	特異度 % (# of Chemical)
EpiOcular*	80 (112)	96 (57)	63 (55)
SkinEthic*	84 (200)	95 (97)	72 (103)
LabCyte**	88 (139)	100 (76)	73 (63)
SkinEthic™ HCE TTT***	87 (151)	94 (94)	75 (57)

Round to the first decimal place when necessary, EpiOcular=EpiOcular™ Eye Irritation Test, SkinEthic=SkinEthic™ Human Corneal Epithelium Eye Irritation Test, LabCyte=LabCyte CORNEA-MODEL24 Eye Irritation Test, #=number

*TG 492²⁾

**Peer Review report, December 2018¹⁾

***TEST GUIDELINE 492B より算出

6. 試験法の適用範囲および留意点

SkinEthic™ HCE TTT は、UN GHS 区分 1、UN GHS 区分 2 への分類および区分に該当しない場合の判定に用いることができるが、UN GHS 区分 2A と 2B を区別することはできない。UN GHS 区分 2A と 2B の区別が必要となる場合は、他の方法やアプローチで対応する必要がある¹²⁾。

本試験法は、プロトコルで定められた濃度に溶解又は均一懸濁した液体の評価が可能である単一物質、混合物、固体、液体、半固体およびワックスに適用することができる。液体は水性と非水性、固体は可溶性と不溶性のいずれも適用可能であるが、水溶性の低い固体 (< 0.014 mg/mL) についてはバリデーション研究において予測性が低かったことから注意が必要である(例：区分 1 の 9 物質のうち 5 つが区分 2 に分類された)。水溶性の低い液体では、20% (w/v) 希釈水溶液で 16 分または 120 分間曝露する際に安定した混合液の状態を保持できない可能性があり、予測性への影響が考えられる。しかしながら、バリデーション試験では水への溶解性が 0.1 g/L 未満の液体 16 物質において予測性に大きな問題はみられなかった。ガスとエアロゾルについてはバリデーション研究で評価されておらず本試験法には適用できない。多成分の物質やその混合物については、技術的には適応可能であるものの、その適応性については限られた情報しか得られていない¹⁾。特定の製剤において本試験が適応できないことを示す根拠がある場合、その製剤は本試験法に適用できない。混合物や試験が困難な化学物質(例：不安定な化学物質)、あるいは適応範囲に明らかに含まれていない物質の試験を検討する場合は、その試験結果が科学的に意味のある結果をもたらすかどうか事前に考慮すべきである。

ホルマザン色素と同じ範囲の光吸収スペクトルを持つ物質や MTT を直接還元する物質は組織生存率の測定に干渉する可能性があり、補正のために適合させた対照物質の使用が必要となる。必要となる対照物質は、被験物質によって生じる干渉の種類および各ホルマザン色素の定量化に使用する手順によって異なる(TG492B の 33-40 項参照)⁶⁾。

7. 結論

SkinEthic™ HCE TTT は、技術移転性、施設内再現性および施設間再現性の高い試験法である。予測性に関しては、OECD 専門家会議が定めた採用基準は満たしていた。

バリデーション報告書の記載に不備は多いが、いずれの施設の記録も十分であり、データも適切であることを PRP が確認しており、バリデーション結果に大きな問題はないと本委員会は考える。

よって、SkinEthic™ HCE TTT 法は、OECD TG492B に準拠して実施した場合、UN GHS 区分 1、UN GHS 区分 2 への分類および区分に該当しない場合の判定を可能とする試験法である。ただし、固体物質の UN GHS 区分 2 の分類は慎重に評価されるべきであると本資料編纂委員会は考える。

8. 参考文献

- 1) OECD Test Guidelines for Chemicals, Available at: <https://www.oecd.org/env/ehs/testing/oecdguidelinesforthetestingofchemicals.htm> (2022 年 6 月現在)
- 2) OECD (2018) Test Guideline 492. Reconstructed human Cornea-like Epithelium (RhCE) test method for identifying chemicals not requiring classification and labelling for eye irritation or serious eye damage.
- 3) Alépée N, Leblanc V, Grandidier MH, Teluob S, Viricel A, Adriaens E, Michaut V.(2021) SkinEthic HCE Time-to-Toxicity on solids: A test method for distinguishing chemicals inducing serious eye damage, eye irritation and not requiring classification and labelling., *Toxicol In Vitro.*,75:105203. doi: 10.1016/j.tiv.2021.105203.
- 4) Alépée N, Leblanc V, Grandidier MH, Teluob S, Tagliati V, Adriaens E, Michaut V.(2020) Development of the SkinEthic HCE Time-to-Toxicity test method for identifying liquid chemicals not requiring classification and labelling and liquids inducing serious eye damage and eye irritation.,*Toxicol In Vitro.*, 269:104960. doi: 10.1016/j.tiv.2020.104960.
- 5) OECD (2022) Test Guideline 492B. Reconstructed human Cornea-like Epithelium (RhCE) test method for Eye Hazard
- 6) Alépée N, Grandidier MH, Teluob S, Amaral F, Caviola E, De Servi B, Martin S, Meloni M, Nardelli L, Padelou C, Tagliati V, Viricel A, Adriaens E, Michaut V. (2022) Validation of the SkinEthic HCE Time-to-Toxicity test method for eye hazard classification of chemicals according to UN GHS.*Toxicol In Vitro.* 80:105319. doi: .1016/j.tiv.2022.105319.
- 7) Barroso, J., Eskes, C., Germolec, D., Kim, B-H., Kojima, H. (2021). Peer Review Panel Report on the scientific validity of the SkinEthic™ Human Corneal Epithelium (HCE) Time-To-Toxicity Test (TTT) on liquids (TTL) and on solids (TTS). January, 12, 2021.
- 8) SkinEthic™ HCE Time-to-Toxicity on Liquids SOP, Version 1. (June 16, 2020). *In Vitro* Eye Irritation Time-To-Toxicity Test Method On Liquids (TTL): SkinEthic™ Human Corneal Epithelium Model (HCE). Available at: [www.episkin.com].
- 9) SkinEthic™ HCE Time-to-Toxicity on Solids SOP, Version 1. (June 16, 2020). *In Vitro* Eye Irritation Time-To-Toxicity Test Method On Solids (TTS): SkinEthic™ Human Corneal Epithelium Model (HCE). Available at: [www.episkin.com].
- 10) OECD (2017) Performance Standards for the Assessment of proposed Similar or Modified *In Vitro* Reconstructed Human Cornea-like Epithelium (RHCE) Test Methods for Eye Hazard
- 11) OECD (2018) PEER REVIEW REPORT OF THE VALIDATION OF THE LABCYTE CORNEA-MODEL²⁴ EYE IRRITATION TEST FOR INCLUSION IN TEST GUIDELINE 492, SERIES ON TESTING AND ASSESMENT No. 282
- 12) OECD (2019) Guidance Document on an Integrated Approach on Testing and Assessment (IATA) for Serious Eye Damage and Eye Irritation Series on Testing & Assessment No. 263 (Second Edition)

Appendix 1

SkinEthic™ HCE TTL (Table 2A) および SkinEthic™ HCE TTS (Table 2B)の 習熟度確認物質リスト

Appendix 1-1 SkinEthic™ HCE TTL リスト

物質名	CASRN	有機官能グループ ¹	物性	生存率1 5分 (%) ^{2,3}	生存率2 16分 (%) ^{2,3}	生存率3 120分 (%) ^{2,3}	VRM 予測性
<i>In Vivo</i> カテゴリー-1							
N,N-Diethylethanolamine	100-37-8	Alcohol. Aliphatic Amine. tertiary	L	2.9±1.7	1.0±1.3	0.6±0.5	Cat 1
Acetic acid (10%)	64-19-7	Carboxylic acid	L	4.2±0.4	25.9±12.7	2.8±0.3	Cat 1
<i>In Vivo</i> カテゴリー-2							
2-Butanone	78-93-3	Ketone	L	21.1±5.4	90.8±13.7	87.6±12.3	Cat 2
Acetone	67-64-1	Ketone	L	6.4±2.6	97.1±3.9	99.1±3.5	Cat 2
Hexadecyltrimethylammonium chloride (2%)	112-02-7	Ammonium salt. Alkylammonium salt	L	58.0±10.7	62.1±9.2	2.4±0.8	Cat 2
<i>In Vivo</i> 非カテゴリー							
1,3-Diisopropylbenzene	99-62-7	Cyclic. Phenyl. Aromatic	L	100.5±8.5	94.5±6.3	96.5±9.3	No Cat
Dodecane	112-40-3	Methyl. Methylene	L	98.3±6.9	103.1±5.6	98.9±5.1	No Cat

Appendix 1-2 SkinEthic™ HCE TTSのリスト

物質名	CASRN	有機官能グループ ¹	物性	生存率1 30分 (%) ^{2,3}	生存率2 120分 (%) ^{2,3}		VRM 予測性
<i>In Vivo</i> カテゴリー-1							
1-Naphthalene acetic acid Na salt	61-31-4	Benzyl. Carboxylic aromatic. Naphthalene	S	1.6±0.7	1.6±0.4		Cat 1
1,2-Benzisothiazol-3(2H)-one	2634-33-5	Benzothiazole/ Benzisothiazole	S	4.0±1.0	4.0±0.4		Cat 1
<i>In Vivo</i> カテゴリー-2							
4-Carboxybenzaldehyde	619-66-9	Aldehyde. Aryl. Carboxylic acid	S	80.1±12.3	5.0±1.4		Cat 2
2-Hydroxy-1,4-naphthoquinone	83-72-7	Enol. Naphthoquinone	S	82.1±8.9	4.5±1.5		Cat 2
Ammonium nitrate	6484-52-2	Amine. Nitrate. Ammonium	S	61.0±4.8	1.8±0.4		Cat 2
<i>In Vivo</i> 非カテゴリー							
Magnesium carbonate	56378-72-4	Magnesium. Carbonate	S	98.3±7.1	93.4±11.7		No Cat
Anthracene	120-12-7	Anthracene. Carboxylic aromatic	S	97.7±6.6	98.1±4.7		No Cat

略号: CASRN = Chemical Abstracts Service Registry Number

VRM = Validated Reference Method, SkinEthic™ HCE TTT

¹ 有機機能グループ OECD QSAR Toolbox 分析に従い分類された (version 3.2; <https://www.oecd.org/chemicalsafety/risk-assessment/oecd-qsar-toolbox.htm>).

² SkinEthic™ HCE TTT バリデーション研究を基に得られた。

³ 液体の場合, the SkinEthic™ HCE TTL は 3 曝露回数を用いる。原体 5 分、20%水溶液 16 分および 120 分である。固体の場合, the SkinEthic™ HCE TTS は 2 曝露回数を用いる。原体 30 分および 120 分である。

添付資料 3

OECD GUIDELINE FOR THE TESTING OF CHEMICALS

Reconstructed Human Cornea-like Epithelium (RHCE) Test Method for Eye Hazard Identification

INTRODUCTION

1. Serious eye damage refers to the production of tissue damage in the eye, or serious physical decay of vision, which is not fully reversible, occurring after exposure of the eye to a test chemical, as defined by the United Nations Globally Harmonized System of Classification and Labelling of Chemicals (UN GHS) (1). Also according to UN GHS, eye irritation refers to the production of changes in the eye, which are fully reversible, occurring after exposure of the eye to a test chemical. Test chemicals inducing serious eye damage are classified as UN GHS Category 1, while those inducing eye irritation are classified as UN GHS Category 2. Test chemicals not classified for eye irritation or serious eye damage are defined as those that do not meet the requirements for classification as UN GHS Category 1 or 2 (2A or 2B) i.e., they are referred to as UN GHS No Category (No Cat).

2. The assessment of serious eye damage/eye irritation has initially involved the use of laboratory animals (OECD Test Guideline (TG) 405; adopted in 1981 and revised in 1987, 2002, 2012 and 2017) (2). *In vitro* or *ex vivo* test methods have been adopted as OECD Test Guidelines (TGs) 437 (4), 438 (5), 460 (6), 491 (7), 492 (8), 494 (9) and 496 (10) to identify either chemicals for serious eye damage potential and/or to identify chemicals not requiring classification for eye hazard potential. The choice of the most appropriate test method to be used should be considered in the context of the OECD Guidance Document on an Integrated Approaches on Testing and Assessment (IATA) for Serious Eye Damage and Eye irritation (3).

3. This TG describes an *in vitro* procedure allowing the identification on its own of chemicals (substances and mixtures) not requiring classification (No Cat), requiring classification for eye irritation (Cat 2) and requiring classification for serious eye damage (Cat 1) according to the UN GHS ocular hazard categories (1).

4. This TG describes a validated test method, namely the SkinEthic™ Human Corneal Epithelium (HCE) Time-to-Toxicity (TTT) test using a commercially available reconstructed Human Cornea-like Epithelium (HCE). The RhCE is designed to closely mimic the histological, morphological, biochemical and physiological properties of the human corneal epithelium. A validation study for assessing the three UN GHS ocular hazard categories has been conducted (11)(12)(13) on this method, referred to in the following text as the Validated Reference Method (VRM). From the validation study and its independent peer review (14) it was concluded that the SkinEthic™ HCE TTT is able to correctly identify chemicals (both substances and mixtures) by discriminating the three UN GHS categories for serious eye damage/eye irritation, i.e. UN GHS Cat. 1, Cat. 2 and No Cat chemicals (1), and the test method was recommended as a full replacement to the *in vivo* Draize acute eye irritation test for classification of chemicals. It is recognized

© OECD, (2022)

You are free to use this material subject to the terms and conditions available at <http://www.oecd.org/termsandconditions/>.

that the use of this test guideline is subject to national and international regulatory considerations and conditions. The Guidance Document No. 263 on IATA should be consulted for further testing with other adequate *in vitro* tests in a weight-of-evidence approach, if deemed necessary (3). Annexes II-V provide a synopsis of the important elements of the test method, as well as flowcharts providing guidance for specific situations.

5. The purpose of this TG is to describe the procedure used to evaluate the eye hazard potential of a test chemical based on its ability to induce cytotoxicity in a RhCE tissue construct, as measured by reduction of a vital dye (MTT [3-(4,5-Dimethylthiazol-2-yl)-2,5-diphenyltetrazolium bromide; Thiazolyl blue tetrazolium bromide; CAS RN 298-93-1]) hereafter designated as the tetrazolium dye (TD) (15) (see paragraph 23). The viability of the RhCE tissue following exposure to a test chemical is determined in comparison to tissues treated with the negative control substance (%viability) for the two or three exposure times, and is then used to predict the eye hazard potential of the test chemical.

6. Definitions are provided in Annex I.

INITIAL CONSIDERATIONS AND LIMITATIONS

7. The SkinEthic™ HCE tissue construct is a three-dimensional model produced using cells from the species of interest which mimics the *in vivo* corneal epithelium(17). The test method directly measures cytotoxicity resulting from penetration of the chemical through the corneal epithelium and production of cell and tissue damage following chemical exposure, which is used to predict the eye hazard identification of a test chemical. Cell damage can occur by several modes of action (see paragraph 16), but cytotoxicity plays an important, if not the primary, mechanistic role in determining the overall serious eye damage/eye irritation response of a chemical, manifested *in vivo* mainly by corneal opacity, iritis, conjunctival redness and/or conjunctival chemosis, regardless of the physicochemical processes underlying tissue damage.

8. A total of 151 chemicals covering a variety of chemical types, chemical classes, chemical structures, as well as molecular weights, LogP and other physical-chemical properties have been tested in the validation study underlying this TG. The validation database covered 134 different organic functional groups (11)(12)(13) and all key *in vivo* drivers of classification (26)(27). The majority of these chemicals represented mono-constituent substances (a total of 151 substances, of which 16 were tested in dilution), but several multi-constituent substances (including surfactants or polymers) were also included in the study. In terms of physical state and according to UN GHS Categories (1), the 151 tested chemicals were distributed as follows: 70 liquids, comprising 21 Cat. 1, 25 Cat. 2 (incl. 16 Cat. 2A and 9 Cat. 2B) and 24 No Cat, and 81 solids, comprising 29 Cat. 1, 19 Cat.2, and 33 No Cat. (11)(12)(13).

9. The SkinEthic™ HCE TTT is not intended to discriminate between UN GHS Category 2A (eye irritation, effects fully reversible within 21 days) and UN GHS Category 2B (mild eye irritation, effects fully reversible within 7 days). This differentiation needs to be addressed by other methods or approaches, if discrimination is deemed necessary (3).

10. This TG is applicable to substances and mixtures, to solids, liquids, semi-solids and waxes. The liquids may be aqueous or non-aqueous, solids may be soluble or insoluble in water. Caution should be used when testing solid chemicals with poor water solubility (< 0.014 mg/mL) as they were frequently underpredicted by SkinEthic™ HCE TTS in the validation study (e.g., 5 out of 9 Cat.1 chemicals were underpredicted to be Cat.2). Gases and aerosols have not been assessed in a validation study. While it is conceivable that these can be tested using RhCE technology, the current TG does not allow testing of gases and aerosols. Limited information is currently available on the applicability of the test method to multi-constituent substances/mixtures (11). The test method is nevertheless technically applicable to the testing of multi-constituent substances and mixtures. In cases where evidence demonstrates the non-

applicability of the TG to specific formulations, the TG should not be used for those formulations. When considering testing of difficult-to-test chemicals (e.g. unstable substances), or test chemicals or mixtures not clearly within the applicability domain described in this Guideline, upfront consideration should be given to whether the results of such testing will yield results that are meaningful scientifically.

11. Test chemicals absorbing light in the same range as formazan dye (FD, naturally or after treatment) and test chemicals able to directly reduce the vital dye MTT (to FD) may interfere with the tissue viability measurements and need the use of adapted controls for corrections. The type of adapted controls that may be required will vary depending on the type of interference produced by the test chemical and the procedure used to quantify each FD (see paragraphs 33-38).

12. The validation study carried out in three laboratories demonstrated that SkinEthic™ HCE TTT is transferable to laboratories considered to be naïve in the conduct of the assays and also is reproducible within- and between laboratories (13)(14). The within-laboratory reproducibility (WLR) for SkinEthic™ was 85-95% for TTL (20 chemicals) and 100% for TTS (20 chemicals). The between-laboratory reproducibility (BLR) was 90% for TTL and 100% for TTS.

13. The SkinEthic™ HCE TTT test can be used to identify chemicals that do not require classification for eye irritation or serious eye damage according to the UN GHS classification system (1). Considering the data obtained (13)(14) (Table 1), the SkinEthic™ HCE TTT test has a balanced accuracy of 74.4% (based on 151 chemicals) with correct predictions of 79% for Cat 1 (based on 50 chemicals), 69% for Cat 2 (based on 44 chemicals) and 75% for No Cat (based on 57 chemicals), when compared to reference *in vivo* rabbit eye test data (OECD TG 405) (2)(16) and classified according to the UN GHS classification system (1). The performance for SkinEthic™ HCE TTL and TTS is different (Annex V).

Table 1. The weighted¹ performance of the SkinEthic™ HCE TTT using the 3x3 matrix showing correct, under- and over- predictions per UN GHS category

UN GHS categories	SkinEthic HCE TTT - Predicted categories (n/N%)		
	Cat 1 (n)	Cat 2 (n)	No Cat (n)
Cat 1 (N=50)	79.2% (39.60)	20.8% (10.40)	0% (0.00)
Cat 2 (N=44)	18.3% (8.06)	69.2% (30.46)	12.5% (5.48)
No Cat (N=57)	1.8% (1.00)	23.3% (13.33)	74.9% (42.67)

14. The term "test chemical" is used in this TG to refer to what is being tested and is not related to the applicability of the RhCE test method to the testing of substances and/or mixtures.

PRINCIPLE OF THE TEST

15. The test chemical is applied topically to a minimum of two three-dimensional RhCE tissue constructs and tissue viability is measured following exposure and a post-soak incubation period. The SkinEthic™ HCE tissues are reconstructed from human immortalized corneal epithelial cells, which have been cultured for several days to form a stratified, highly differentiated squamous epithelium morphologically similar to that found in the human cornea. The SkinEthic™ HCE tissue construct consists

¹ In the weighted calculation, each chemical has an equal weight of 1 in the performance regardless of the number of test runs (i.e., in 3 overall test runs, if chemical A is tested Cat.1 two times and Cat.2 one time, a fractional weight of 0.66 (2/3) and 0.33 (1/3) is assigned to the number of predictions (n) for Cat.1 and Cat.2, respectively (11).

of at least 4 viable layers of cells including columnar basal cells, transitional wing cells and superficial squamous cells similar to that of the normal human corneal epithelium (17)(18).

16. Chemical-induced serious eye damage/eye irritation, manifested *in vivo* mainly by corneal opacity, iritis, conjunctival redness and/or conjunctival chemosis, is the result of a cascade of events beginning with penetration of the chemical through the cornea and/or conjunctiva and production of damage to the cells. Cell damage can occur by several modes of action, including: cell membrane lysis (e.g., by surfactants, organic solvents); coagulation of macromolecules (particularly proteins) (e.g., by surfactants, organic solvents, alkalis and acids); saponification of lipids (e.g., by alkalis); and alkylation or other covalent interactions with macromolecules (e.g., by bleaches, peroxides and alkylators) (19)(20)(21). However, it has been shown that cytotoxicity plays an important, if not the primary, mechanistic role in determining the overall serious eye damage/eye irritation response of a chemical regardless of the physicochemical processes underlying tissue damage (22). Moreover, the serious eye damage/eye irritation potential of a chemical is principally determined by the extent of initial injury (20), which correlates with the extent of cell death (22) and with the extent of the subsequent responses and eventual outcomes (23)(24). Thus, slight irritants generally only affect the superficial conjunctival and corneal epithelium, the mild and moderate irritants damage principally the epithelium and superficial stroma and the severe irritants damage the epithelium, deep stroma and at times the endothelium (22)(25). The measurement of viability of the RhCE tissue construct after topical exposure to a test chemical to identify chemicals not requiring classification for serious eye damage/eye irritancy (UN GHS No Category) is based on the assumption that all chemicals inducing serious eye damage or eye irritation will induce cytotoxicity in the corneal epithelium and/or conjunctiva.

17. RhCE tissue viability is classically measured by enzymatic conversion of MTT for the SkinEthic™ HCE TTT by the viable cells of the tissue into a formazan dye (FD), which is quantitatively measured after extraction from tissues (15).

18. The SkinEthic™ HCE TTT Test is based on two protocols, one for liquids (SkinEthic™ HCE TTL) and one for solids (SkinEthic™ HCE TTS). The SkinEthic™ HCE TTT uses three exposure times for TTL and two exposure times for TTS (see paragraphs 29-30). The SkinEthic™ HCE TTL and TTS protocols make use of different prediction models. For SkinEthic™ HCE TTL, a chemical that results in a mean viability below or equal to 50% within all time treatments will be classified as a Cat 1, and one that results in a mean viability strictly above 50% as a No Cat. Any other combination of mean viability values will lead to a Cat 2 classification (see paragraph 43). For SkinEthic™ HCE TTS, a chemical that results in a mean viability below or equal to 40% after 30 minutes. exposure and below or equal to 60% after 120 minutes. exposure will be classified as a Cat 1. A mean viability strictly above these cut-offs within the two time treatments will classify the chemical as a No Cat. Any other combination of mean viability values will classify the chemical as a Cat 2 (see paragraph 43).

DEMONSTRATION OF PROFICIENCY

19. Prior to routine use of the test method for regulatory purposes, laboratories should demonstrate technical proficiency by correctly predicting the fourteen proficiency chemicals listed in Table 2. These chemicals were selected from the chemicals used in the validation studies of the SkinEthic™ HCE TTT (13)(14). The selection includes, to the extent possible, chemicals that: (i) cover different physical states; (ii) cover the range of *in vivo* serious eye damage/eye irritation responses based on high quality results obtained in the reference *in vivo* rabbit eye test (OECD TG 405) (2)(16) and the UN GHS classification system (i.e., Categories 1, 2, or No Category) (1); (iii) cover the key *in vivo* drivers of classification (26)(27); (iv) are representative of the chemical classes used in the validation study (13); (v) cover a good and wide representation of organic functional groups (11)(12)(13); (vi) have chemical structures that are well-defined (11)(12)(13); (vii) produced reproducible results in the RhCE test method during the validation study; (viii)

were correctly predicted by the RhCE test method during the validation study; (ix) cover the range of *in vitro* responses based on high quality of the test method data; (x) are commercially available; and (xi) are not associated with prohibitive acquisition and/or disposal costs. In situations where a listed chemical is unavailable or cannot be used for other justified reasons, another chemical fulfilling the criteria described above, e.g. from the chemicals used ² the validation of the SkinEthic™ HCE TTT, could be used. Such deviations should however be justified.

Table 2: List of proficiency chemicals for SkinEthic™ HCE TTL (Table 2A) and for SkinEthic™ HCE TTS (Table 2B)

Table 2A. List of Proficiency Chemicals for SkinEthic™ HCE TTL

ChemicalName	CASRN	OrganicFunctionalGroup ¹	Physicalstate	Viability 1 5 min. (%) ^{2,3}	Viability 2 16 min. (%) ^{2,3}	Viability 3 120 min. (%) ^{2,3}	VRM Prediction
<i>In Vivo</i>Category1							
N.N-Diethylethanolamine	100-37-8	Alcohol. Aliphatic Amine. tertiary	L	2.9±1.7	1.0±1.3	0.6±0.5	Cat 1
Acetic acid (10%)	64-19-7	Carboxylic acid	L	4.2±0.4	25.9±12.7	2.8±0.3	Cat 1
<i>In Vivo</i>Category2							
2-Butanone	78-93-3	Ketone	L	21.1±5.4	90.8±13.7	87.6±12.3	Cat 2
Acetone	67-64-1	Ketone	L	6.4±2.6	97.1±3.9	99.1±3.5	Cat 2
Hexadecyltrimethylammonium chloride (2%)	112-02-7	Ammonium salt. Alkylammonium salt	L	58.0±10.7	62.1±9.2	2.4±0.8	Cat 2
<i>In Vivo</i>NoCategory							
1,3-Diisopropylbenzene	99-62-7	Cyclic. Phenyl. Aromatic	L	100.5±8.5	94.5±6.3	96.5±9.3	No Cat
Dodecane	112-40-3	Methyl. Methylene	L	98.3±6.9	103.1±5.6	98.9±5.1	No Cat

Table 2B. List of Proficiency Chemicals for SkinEthic™ HCE TTS

ChemicalName	CASRN	OrganicFunctionalGroup ¹	Physicalstate	Viability 1 30 min. (%) ^{2,3}	Viability 2 120 min. (%) ^{2,3}	VRM Prediction
<i>In Vivo</i>Category1						
1-Naphthalene acetic acid Na salt	61-31-4	Benzyl. Carboxylic aromatic. Naphthalene	S	1.6±0.7	1.6±0.4	Cat 1
1,2-Benzisothiazol-3(2H)-one	2634-33-5	Benzothiazole/ Benzoisothiazole	S	4.0±1.0	4.0±0.4	Cat 1
<i>In Vivo</i>Category2						
4-Carboxybenzaldehyde	619-66-9	Aldehyde. Aryl. Carboxylic acid	S	80.1±12.3	5.0±1.4	Cat 2
2-Hydroxy-1,4-naphthoquinone	83-72-7	Enol. Naphthoquinone	S	82.1±8.9	4.5±1.5	Cat 2
Ammonium nitrate	6484-52-2	Amine. Nitrate. Ammonium	S	61.0±4.8	1.8±0.4	- Cat 2
<i>In Vivo</i>NoCategory						
Magnesium carbonate	56378-72-4	Magnesium. Carbonate	S	98.3±7.1	93.4±11.7	No Cat
Anthracene	120-12-7	Anthracene. Carboxylic aromatic	S	97.7±6.6	98.1±4.7	No Cat

1. 3-Diisopropylbenzene

2. 2-Benzisothiazol-3(2H)-one

Abbreviations: CASRN = Chemical Abstracts Service Registry Number; VRM = Validated Reference Method, SkinEthic™ HCE TTT.

¹Organic functional group assigned according to an OECD QSAR Toolbox analysis (version 3.2;<https://www.oecd.org/chemicalsafety/risk-assessment/oecd-qsar-toolbox.htm>).

²Based on results obtained with SkinEthic™ HCE TTT in the validation study (13)(14).

³For liquids, the SkinEthic™ HCE TTL uses three exposure-times : Test chemicals are applied neat for 5 minutes (viability 1), and diluted at 20% (w/v) for 16 minutes (viability 2) and 120 minutes (viability 3). For solids, the SkinEthic™ HCE TTS uses two exposure-times: Test chemicals are applied neat for 30 minutes (viability 1) and 120 minutes (viability 2) (see paragraph 29).

20. As part of the proficiency testing, it is recommended that users verify the barrier properties of the tissues after receipt as specified by the RhCE tissue construct producer (see paragraphs 24, 25 and 27). This is particularly important if tissues are shipped over long distance / time periods. Once a test method has been successfully established and proficiency in its use has been acquired and demonstrated, such verification will not be necessary on a routine basis. However, when using a test method routinely, it is recommended to continue to assess the barrier properties at regular intervals.

PROCEDURE

21. The test method currently covered by this TG is the scientifically validated SkinEthic™ HCE TTT Test (14). The Standard Operating Procedures (SOP) for both SkinEthic™ HCE TTL and TTS are available and should be employed when implementing and using the test method in a laboratory (28)(29). The following paragraphs and Annex II describe the main components and procedures of the RhCE test method.

RhCE TEST METHOD Components

GENERAL CONDITIONS

22. Relevant human-derived cells should be used to reconstruct the cornea-like epithelium three-dimensional tissue, which should be composed of progressively stratified but not cornified cells. The RhCE tissue construct is prepared in inserts with a porous synthetic membrane through which nutrients can pass to the cells. Multiple layers of viable, non-keratinized epithelial cells should be present in the reconstructed cornea-like epithelium. The RhCE tissue construct should have the epithelial surface in direct contact with air so as to allow for direct topical exposure of test chemicals in a fashion similar to how the corneal epithelium would be exposed *in vivo*. The RhCE tissue construct should form a functional barrier with sufficient robustness to resist rapid penetration of cytotoxic benchmark substances, e.g., sodium dodecyl sulphate (SDS). For example, the barrier function should be demonstrated and may be assessed by determination of the concentration at which a benchmark substance reduces the viability of the tissues by 50% (IC₅₀) following a fixed exposure time e.g., 30 minutes-treatment with 50 µL SDS (see paragraph 27). The containment properties of the RhCE tissue construct should prevent the passage of test chemical around the edge of the viable tissue, which could lead to poor modelling of corneal exposure. The human-derived cells used to establish the RhCE tissue construct should be free of contamination by bacteria, viruses, mycoplasma, and fungi. The sterility of the tissue construct should be checked by the supplier for absence of contamination by fungi and bacteria.

FUNCTIONAL CONDITIONS

Viability

23. The assay used for quantifying tissue viability is the tetrazolium dye (MTT) assay (15). Viable cells of the RhCE tissue construct reduce the vital dye MTT into a blue MTT formazan precipitate, which is then extracted from the tissue using isopropanol (or a similar solvent). The extracted formazan dye may be quantified using either a standard absorbance (Optical Density (OD)) measurement or an HPLC/UPLC-spectrophotometry procedure (28)(29). The OD of the blank solution alone (which is the extraction solvent for MTT assay) should be sufficiently small, i.e., $OD < 0.1$. Users of the RhCE tissue construct should ensure that each batch of the RhCE tissue construct used meets defined criteria for the negative control. Acceptability ranges for the negative control OD values for the SkinEthic™ HCE TTT are given in Table 3. For HPLC/UPLC-spectrophotometry, the negative control OD ranges provided in Table 3 should be used as acceptance criterion for the negative control. It should be documented in the test report that the tissues treated with the negative control substance are stable in culture (provide similar tissue viability measurements) for the duration of the test exposure period. A similar procedure should be followed by the tissue producer as part of the quality control tissue batch release, but in this case different acceptance criteria than those specified in Table 3 may apply. An acceptability range (upper and lower limit) for the negative control OD values (in the QC test method conditions) should be established by the RhCE tissue construct developer/supplier.

Table 3. Acceptability ranges for negative control OD values (for the test method users)

Test Method	Lower acceptance limit	Upper acceptance limit
SkinEthic™ HCE TTT(HCE/S) - VRM (for both TTL and TTS protocols)	> 1.0	≤ 2.5

Barrier function

24. The RhCE tissue construct should be sufficiently thick and robust to resist the rapid penetration of cytotoxic benchmark substances, as estimated e.g. by IC_{50} (SDS) (Table 4). The barrier function of each batch of the RhCE tissue construct used should be demonstrated by the RhCE tissue construct developer/vendor upon supply of the tissues to the end user (see paragraph 27).

Morphology

25. Histological examination of the RhCE tissue construct should demonstrate human cornea-like epithelium structure (including at least four layers of viable epithelial cells). For the VRM, appropriate morphology has been established by the supplier and therefore does not need to be demonstrated again by a test method user for each tissue batch used.

Reproducibility

26. The results of the positive and negative controls of the test method should demonstrate reproducibility over time.

Quality control (QC)

27. The RhCE tissue construct should only be used if the developer/supplier demonstrates that each batch of the RhCE tissue construct used meets defined production release criteria, the most relevant among which are those for viability (paragraph 23) and barrier function. An acceptability range (upper and lower limits) for the barrier functions as measured by the IC₅₀ should be established by the RhCE tissue construct supplier. The IC₅₀ acceptability range used as QC batch release criterion by the supplier of the RhCE tissue constructs used in the test method is given in Table 4. Data demonstrating compliance with all production release criteria should be provided by the RhCE tissue construct supplier to the test method users so that they are able to include this information in the test report. Only results produced with tissues fulfilling all of these production release criteria can be accepted for reliable prediction of chemicals for eye hazard identification in accordance with UN GHS.

Table 4. QC batch release criteria of the SkinEthic™ HCE tissue construct

Test Method	Lower acceptance limit	Upper acceptance limit
SkinEthic™ HCE TTT(HCE/S) - VRM (30 minutes treatment with 50 µL SDS)	IC ₅₀ = 1.0 mg/mL	IC ₅₀ = 3.2 mg/mL

Application of the Test Chemical and Control Substances

28. At least two tissue replicates for each time-treatment point should be used for each test chemical and each control substance in each run. Two different treatment protocols are used, one for liquid test chemicals and one for solid test chemicals (28)(29). For the SkinEthic™ HCE TTT, the tissue construct surface should be moistened with calcium and magnesium-free Dulbecco's Phosphate Buffered Saline (Ca²⁺/Mg²⁺-free DPBS) for liquids' protocol before application of test chemicals or with water for solids' protocol, to mimic the wet conditions of the human eye. The treatment of the tissues is initiated with exposure to the test chemical(s) and control substances. For any treatment protocols, a sufficient amount of test chemical or control substance should be applied to uniformly cover the epithelial surface.

29. Test chemicals that can be pipetted at 37°C or lower (using a positive displacement pipette, if needed) are treated as liquids in the test method, otherwise they should be treated as solids (see paragraph 30). Specific procedures for liquids that do not form a solution or a stable suspension are included in the SOP (28). In the test method, liquid test chemicals are evenly spread over the tissue surface (i.e. 160±2 µL/cm² application) (Annex II). SkinEthic™ HCE tissues are topically exposed for three different treatment times: Test chemicals are applied neat for 5 minutes, and diluted at 20% (w/v) in distilled water for 16 minutes and 120 minutes at the pre-defined conditions of the method (28). At the end of the exposure period, the liquid test chemical and the control substances should be carefully removed from the tissue surface by extensive rinsing with Ca²⁺/Mg²⁺-free DPBS at room temperature. This rinsing step should be followed by a 10 minute post-exposure immersion in fresh medium at room temperature (to remove any test chemical absorbed into the tissue) prior to performing the MTT assay (Annex II)(28).

30. Test chemicals that cannot be pipetted (for example, highly viscous liquids) are treated as solids in the test method. The amount of test chemical applied should be sufficient to cover the entire surface of the tissue, i.e. a minimum of 160±2 mg/cm² application should be used (Annex II). Whenever possible, solids should be tested as a fine powder. SkinEthic™ HCE tissues are topically exposed to test chemicals for 30 and 120 minutes at standard culture conditions (see Annex II)(29). At the end of the exposure period, the solid test chemical and the control substances should be carefully removed from the tissue surface by extensive rinsing with Ca²⁺/Mg²⁺-free DPBS at room temperature. This rinsing step should be followed by

a 30 minute post-exposure immersion in fresh medium at room temperature (to remove any test chemical absorbed into the tissue), prior to performing the TD assay (Annex II)(29).

31. Concurrent negative and positive controls should be included in each run to demonstrate that the viability (determined with the negative control) and the sensitivity (determined with the positive control) of the tissues are within acceptance ranges defined based on historical data. The concurrent negative control also provides the baseline (100% tissue viability) to calculate the relative percent viability of the tissues treated with the test chemical (%Viability_{test}). The recommended positive control substance to be used with the SkinEthic™ HCE TTT is neat methyl acetate (CAS RN. 79-20-9, commercially available) for liquids' protocol, and lactic acid at 1% (W/V) in water (CAS RN. 50-21-5, commercially available) for solids' protocol. The recommended negative control substance to be used is Ca²⁺/Mg²⁺-free DPBS for both liquids and solids protocols. These were the control substances used in the validation studies of the SkinEthic™ HCE TTT and are those for which most historical data exist. The use of suitable alternative positive or negative control substances should be scientifically and adequately justified. Negative and positive controls should be tested with the same protocol(s) as the one(s) used for the test chemicals included in the run (i.e. for liquids and/or solids). This application should be followed by the treatment exposure, rinsing, and a post-exposure immersion, as described for controls run concurrently to liquid test chemicals (see paragraph 29) or for controls run concurrently to solid test chemicals (see paragraph 30), prior to performing the MTT assay (see paragraph 32) (26)(27). One single set of negative and positive controls is sufficient for all test chemicals of the same physical state (liquids or solids) included in the same run.

Tissue Viability Measurements

32. The MTT assay is a standardised quantitative method (15) that should be used to measure tissue viability under this TG. It is compatible with use in a three-dimensional tissue construct. The MTT assay is performed immediately following the post-soak procedure. The RhCE tissue construct sample is placed in 0.3 mL of MTT solution at 1 mg/mL for 180±15 minutes at standard culture conditions. The vital dye MTT is reduced into a blue MTT formazan precipitate by the viable cells of the RhCE tissue construct. The precipitated blue MTT formazan product is then extracted from the tissue using an appropriate volume of isopropanol (or a similar solvent) (28)(29). Tissues tested with liquid test chemicals should be extracted from both the top and the bottom of the tissues, while tissues tested with solid test chemicals or coloured liquids should be extracted from the bottom of the tissue only (to minimise any potential contamination of the isopropanol extraction solution with any test chemical that may have remained on the tissue). Tissues tested with liquid test chemicals that are not readily washed off may also be extracted from the bottom of the tissue only. The concurrently tested negative and positive control substances should be treated similarly to the tested chemical. The extracted MTT formazan may be quantified either by a standard absorbance (OD) measurement at 570 nm using a filter band pass of maximum ± 30 nm or by using an HPLC/UPLC-spectrophotometry procedure (see paragraph 40).

33. Optical properties of the test chemical or its chemical action on MTT may interfere with the measurement of formazan dye (FD) leading to a false estimate of tissue viability, i.e., under-prediction of eye irritation. Test chemicals may interfere with the measurement of FD by direct reduction of the MTT into coloured FD (blue MTT formazan) and/or by colour interference if the test chemical absorbs, naturally or due to treatment procedures, in the same OD range as FD (i.e., MTT formazan: around 570 nm). Potential of chemicals to directly reduce TD and/or interfere with colour (only necessary for coloured test chemicals) should be checked before testing. In case of FD interference, additional controls should be used to correct for potential interference from such test chemicals (see paragraphs 34-36 and Annexes III). This is especially important when a specific test chemical is not completely removed from the RhCE tissue construct by rinsing or when it penetrates the cornea-like epithelium and is therefore present in the RhCE tissue constructs when the MTT assay is performed. For test chemicals absorbing light in the same range

as FD (naturally or after treatment), which are not compatible with the standard absorbance (OD) measurement of FD due to strong interference, i.e., strong absorption at 570 ± 30 nm (with MTT formazan), an HPLC/UPLC-spectrophotometry procedure to measure FD may be employed (see paragraphs 39 and 40) (28). A detailed description of how to detect and correct for direct MTT reduction and interferences by colouring agents is available in the test method SOPs (28)(29). Illustrative flowcharts providing guidance on how to identify and handle direct MTT-reducers and/or colour interfering chemicals are also provided in Annexes III.

34. To identify potential interference by test chemicals absorbing light in the same range as FD (naturally or after treatment) and decide on the need for additional controls, the test chemical is added to water and incubated for an appropriate time at room temperature (Annex II) (28)(29). If a coloured solution is obtained when mixing the test chemical with water (see Annex III), the test chemical is presumed to interfere with the standard absorbance (OD) measurement of FD; in such a case, further colorant controls should be performed or, alternatively, an HPLC/UPLC-spectrophotometry procedure should be used, which do not require these controls (see paragraphs 39 and 40 and Annexes III and IV). When performing the standard absorbance (OD) measurement, each interfering test chemical should be applied on at least two viable tissue replicates, which undergo the entire testing procedure but are incubated with medium instead of MTT solution during the MTT incubation step, to generate a non-specific colour in living tissues (NSC_{living}) control (28)(29). The NSC_{living} control needs to be performed concurrently to the testing of the coloured test chemical and, in case of multiple testing, an independent NSC_{living} control needs to be conducted with each test performed (in each run) due to the inherent biological variability of living tissues. True tissue viability is calculated as: the percent tissue viability obtained with living tissues exposed to the interfering test chemical and incubated with the MTT (%Viability_{test}) minus the percent non-specific colour obtained with living tissues exposed to the interfering test chemical and incubated with medium without MTT, run concurrently to the test being corrected (%NSC_{living}), i.e., True tissue viability = [%Viability_{test}] - [%NSC_{living}].

35. To identify direct MTT reducers, each test chemical should be added to freshly prepared MTT solution. An appropriate amount of test chemical is added to a MTT solution and the mixture is incubated for approximately 3 hours at standard culture conditions (see Annex III) (28)(29). If the MTT mixture containing the test chemical (or suspension for insoluble test chemicals) turns blue/purple, the test chemical is presumed to directly reduce the MTT and a further functional check on non-viable RhCE tissue constructs should be performed, independently of using the standard absorbance (OD) measurement or an HPLC/UPLC- spectrophotometry procedure. This additional functional check employs killed tissues that possess only residual metabolic activity but absorb and retain the test chemical in a similar way as viable tissues. Killed tissues of SkinEthic™ HCE TTT are prepared by prolonged incubation (e.g., at least 24 ± 1 hours) in water followed by storage at low temperature ("water-killed"). Each MTT reducing test chemical is applied on at least two killed tissue replicates, which undergo the entire testing procedure, to generate a non-specific MTT reduction (NSMTT) control (28)(29). A single NSMTT control is sufficient per test chemical regardless of the number of independent tests/runs performed. True tissue viability is calculated as: the percent tissue viability obtained with living tissues exposed to the MTT reducer (%Viability_{test}) minus the percent non-specific TD reduction obtained with the killed tissues exposed to the same reducer, calculated relative to the negative control run concurrently to the test being corrected (%NSMTT), i.e. true tissue viability = [%Viability_{test}] - [%NSMTT].

36. Test chemicals that are identified as producing both colour interference (see paragraph 34) and direct MTT reduction (see paragraph 35) will also require a third set of controls when performing the standard absorbance (OD) measurement, apart from the NSMTT and NSC_{living} controls described in the previous paragraphs. This is usually the case with darkly coloured test chemicals absorbing light in the range of 570 ± 30 nm for MTT formazan (e.g., blue, purple, black) because their intrinsic colour impedes the assessment of their capacity to directly reduce MTT as described in paragraph 35. This forces the use of NSMTT controls, by default, together with the NSC_{living} controls. Test chemicals for which both NSMTT

and NSC_{living} controls are performed may be absorbed and retained by both living and killed tissues. Therefore, in this case, the NSMTT control may not only correct for potential direct MTT reduction by the test chemical, but also for colour interference arising from the absorption and retention of the test chemical by killed tissues. This could lead to double correction for colour interference since the NSC_{living} control already corrects for colour interference arising from the absorption and retention of the test chemical by living tissues. To avoid a possible double correction for colour interference, a third control for non-specific colour in killed tissues (NSC_{killed}) needs to be performed (Annex III) (28)(29). In this additional control, the test chemical is applied on at least two killed tissue replicates, which undergo the entire testing procedure but are incubated with medium instead of MTT solution during the MTT incubation step. A single NSC_{killed} control is sufficient per test chemical regardless of the number of independent tests/runs performed, but should be performed concurrently to the NSMTT control and with the same tissue batch. True tissue viability is calculated as: the percent tissue viability obtained with living tissues exposed to the test chemical (%Viability_{test}) minus %NSMTT minus %NSC_{living} plus the percent non-specific colour obtained with killed tissues exposed to the interfering test chemical and incubated with medium without MTT, calculated relative to the negative control ran concurrently to the test being corrected (%NSC_{killed}), i.e., True tissue viability = [%Viability_{test}] - [%NSMTT] - [%NSC_{living}] + [%NSC_{killed}].

37. It is important to note that non-specific MTT reduction and non-specific colour interferences may increase the OD (when performing standard absorbance measurements) of the sample above the linearity range of the spectrophotometer and that non-specific MTT reduction can also increase the FD peak area (when performing HPLC/UPLC- spectrophotometry measurements) of the sample above the linearity range of the spectrophotometer. On this basis, when using RhCEs, it is important for each laboratory to determine the OD/peak area linearity range of their spectrophotometer with MTT formazan (CAS RN. 57360-69-7) which is commercially available.

38. The standard absorbance (OD) measurement using a spectrophotometer is appropriate to assess direct MTT-reducers and colour interfering test chemicals, when the observed interference with the measurement of FD is not strong (i.e., the ODs of the samples obtained with the test chemical without any correction for direct MTT reduction and/or colour interference are within the linear range of the spectrophotometer). Nevertheless, results for test chemicals producing %NSMTT and/or %NSC_{living} ≥ 50% for liquids protocol or 60% for solids protocol of the negative control should be taken with caution. Standard absorbance (OD) can however not be measured when the interference with the measurement of FD is strong (i.e., leading to uncorrected ODs of the test samples falling outside of the linear range of the spectrophotometer). Coloured test chemicals or test chemicals that become coloured in contact with water that interfere strongly with the standard absorbance (OD) measurement of FD may still be assessed using HPLC/UPLC-spectrophotometry (Annex III). This is because the HPLC/UPLC system allows for the separation of the FD from the chemical before its quantification. For this reason, NSC_{living} or NSC_{killed} controls are never required when using HPLC/UPLC-spectrophotometry, independently of the chemical being tested. NSMTT control should nevertheless be used if the test chemical is suspected to directly reduce MTT (following the procedure described in paragraph 35). NSMTT control should also be used with test chemicals having a colour (intrinsic or appearing when in water) that impedes the assessment of their capacity to directly reduce MTT as described in paragraph 35. When using HPLC/UPLC-spectrophotometry to measure FD, the percent tissue viability is calculated as percent FD peak area obtained with living tissues exposed to the test chemical relative to the FD peak obtained with the concurrent negative control. For test chemicals able to directly reduce MTT, true tissue viability is calculated as: %Viability_{test} minus %NSMTT as described in the last sentence of paragraph 35.

39. Finally, it should be noted that direct MTT-reducers that are also colour interfering, which are retained in the tissues after treatment and reduce MTT so strongly that they lead to ODs (using standard OD measurement) or peak areas (using UPLC/HPLC- spectrophotometry) of the tested samples that fall outside of the linearity range of the spectrophotometer cannot be assessed with the RhCE test method, although these are expected to occur in only very rare situations.

40. HPLC/UPLC-spectrophotometry may be used with all types of test chemicals (coloured, non-coloured, MTT-reducers and non-MTT reducers) for measurement of FD. Due to the diversity of HPLC/UPLC-spectrophotometry systems, it is not feasible for each user to establish the exact same system conditions. As such, qualification of the HPLC/UPLC-spectrophotometry system should be demonstrated before its use to quantify MTT from samples by meeting the acceptance criteria for a set of standard qualification parameters based on those described in the U.S. Food and Drug Administration guidance for industry on bioanalytical method validation (30)(31). These key parameters and their acceptance criteria are shown in Annex IV. Once the acceptance criteria defined in Annex IV have been met, the HPLC/UPLC-spectrophotometry system is considered qualified and ready to measure FD under the experimental conditions described in this TG.

Acceptance Criteria

41. For each run using RhCE tissue batches that met the quality control (see paragraph 27), tissues treated with the negative control substance should exhibit OD reflecting the quality of the tissues that followed shipment, receipt steps and all protocol processes and should not be outside the historically established boundaries described in Table 3(see paragraph23). For the SkinEthic™ HCE TTL protocol, tissues treated with the positive control substance, i.e., methyl acetate, should show a mean tissue viability $\leq 50\%$ at the time exposure of 5 minutes and $> 50\%$ at both times 16 and 120 minutes, relative to the negative control. For the SkinEthic™ HCE TTS protocol, tissues treated with the positive control substance, i.e., lactic acid at 1% (W/V) in water, should show a mean tissue viability $> 40\%$ relative to the negative control at the time exposure of 30 minutes and $20\% < \%Viability \leq 70\%$ after 120 minutes treatment, thus reflecting the ability of the tissues to respond to an irritant test chemical under the conditions of the test method (28)(29).

42. The variability between tissue replicates of test chemicals and control substances should fall within the accepted limits (i.e., the difference of viability between two tissue replicates should be less than 20% or the standard deviation (SD) between three tissue replicates should not exceed 18%). If either the negative control or positive control included in a run is outside of the accepted ranges, the run is considered "non-qualified" and should be repeated. If the variability between tissue replicates of a test chemical is outside of the accepted range, the test must be considered "non-qualified" and the test chemical should be re-tested.

Interpretation of Results and Prediction Model

43. The OD values/peak areas obtained with the replicate samples for each test chemical should be used to calculate the mean percent tissue viability (mean between tissue replicates) normalised to the negative control, which is set at 100%. The percentage tissue viability cut-off values used to identify and classify test chemicals for ocular hazard categories according to UN GHS classification; i.e. No Category (not classified), Category 2 (eye irritation) and Category 1 (serious eye damage) are given in Table 5. Results should thus be interpreted as follows:

- The test chemical is identified as a chemical not requiring classification and labelling according to UN GHS (No Category) if the mean percent tissue viability after exposure and a post-soak incubation is more than ($>$) the established percentage tissue viability cut-off value within all-time treatments.
- The test chemical is identified as a chemical inducing serious eye damage according to UN GHS (Category 1) if the mean percent tissue viability after exposure and a post-soak incubation is less than or equal to (\leq) the established percentage tissue viability cut-off value within all-time treatments.

- The test chemical is identified as a chemical inducing eye irritation according to UN GHS (Category 2) if the combination of mean percent tissue viability within all-time treatment fall outside of the criteria established to identify the test chemical as a No Cat or a Cat 1, as shown in Table 5.

Table 5. Prediction Model according to UN GHS classification

	No Category	Category 2	Category 1
SkinEthic HCE TTL (for the liquids protocol)	Mean tissue viability > 50% within all-time treatments	Any other combination of values ¹	Mean tissue viability ≤ 50% within all-time treatments
SkinEthic HCE TTS _(for the solids protocol)	Mean tissue viability > 40% after 30 minutes and > 60% after 120 minutes	Any other combination of values ¹	Mean tissue viability ≤ 40% after 30 minutes and ≤ 60% after 120 minutes

¹Any combination of values other than those defined for No Cat or Cat 1.

44. A single test composed of at least two tissue replicates should be sufficient for a test chemical when the result is unequivocal. However, in cases of borderline results, such as non-concordant replicate measurements and/or mean percent tissue viability equal to established percentage tissue viability cut-off value $\pm 5\%$ within the time treatments, a second test should be considered, as well as a third one in case of discordant results between the first two tests.

DATA AND REPORTING

Data

45. Data from individual replicate tissues in a run (e.g., OD values/FD peak areas and calculated percent tissue viability data for the test chemical and controls, and the final RhCE test method prediction) should be reported in tabular form for each test chemical, including data from repeat tests, as appropriate. In addition, mean percent tissue viability and difference of viability between two tissue replicates (if $n=2$ replicate tissues) or SD (if $n \geq 3$ replicate tissues) for each individual test chemical and control should be reported. Any observed interferences of a test chemical with the measurement of FD through direct MTT reduction and/or coloured interference should be reported for each tested chemical.

Test report

46. The test report should include the following information:

Test chemical*Mono-constituent substance*

- Chemical identification, such as IUPAC or CAS name(s), CAS registry number(s), SMILES or InChI code, structural formula, and/or other identifiers;
- Physical state, volatility, pH, LogP, molecular weight, chemical class, and additional relevant physicochemical properties relevant to the conduct of the study, to the extent available;
- Purity, chemical identity of impurities as appropriate and practically feasible, etc.;
- Treatment prior to testing, if applicable (e.g., warming, grinding);
- Storage conditions and stability to the extent available.

Multi-constituent substance, UVCB and mixture

- Characterisation as far as possible by e.g., chemical identity (see above), purity, quantitative occurrence and relevant physicochemical properties (see above) of the constituents, to the extent possible
- Physical state and additional relevant physicochemical properties relevant to the conduct of the study, to the extent possible;
- Purity, chemical identity of impurities as appropriate and practically feasible, etc.;
- Treatment prior to testing, if applicable (e.g., warming, grinding);
- Storage conditions and stability to the extent possible.

Positive and Negative Control Substances

- Chemical identification, such as IUPAC or CAS name(s), CAS registry number(s), SMILES or InChI code, structural formula, and/or other identifiers;
- Physical state, volatility, molecular weight, chemical class, and additional relevant physicochemical properties relevant to the conduct of the study, to the extent available;
- Purity, chemical identity of impurities as appropriate and practically feasible, etc.;
- Treatment prior to testing, if applicable (e.g., warming, grinding);
- Storage conditions and stability to the extent available;
- Justification for the use of a different negative control than those referenced in Annex II, if applicable;
- Justification for the use of a different positive control than those referenced in Annex II, if applicable;
- Reference to historical positive and negative control results demonstrating suitable run acceptance criteria.

Information Concerning the Sponsor and the Test Facility

- Name and address of the sponsor, test facility and study director.

*RhCE Tissue Construct and Protocol Used (providing rationale for the choices)**Test Method Conditions*

- RhCE tissue construct used, including batch number;

- Wavelength and band pass (if applicable) used for quantifying FD, and linearity range of measuring device (e.g., spectrophotometer);
- Description of the method used to quantify FD
- Description of the HPLC/UPLC-spectrophotometry system used, if applicable; LLOQ, ULOQ, and results of calibration curves and QCs using the same type of fitting and weighting as in the validation of the HPLC/UPLC should be included in the reporting.
- Complete supporting information for the specific RhCE tissue construct used including its performance. This should include, but is not limited to:
 - i. Viability quality control (supplier)
 - ii. Viability under test method conditions (user);
 - iii. Barrier function quality control (supplier);
 - iv. Morphology, if available (supplier);
 - v. Other quality controls (QC) of the RhCE tissue construct, if available (supplier);
- Reference to historical data of the RhCE tissue construct. This should include, but is not limited to: Acceptability of the QC data with reference to historical batch data;
- Statement that the testing facility has demonstrated proficiency in the use of the test method before routine use by testing of the proficiency chemicals;

Run and Test Acceptance Criteria

- Positive and negative control means and acceptance ranges based on historical data;
- Acceptable variability between tissue replicates for positive and negative controls;
- Acceptable variability between tissue replicates for the test chemical;

Test Procedure:

- Details of the test procedure used (e.g., version of the SOP);
- Doses of test chemical and control substances used;
- Duration and temperature of exposure, post-exposure immersion and post-exposure incubation periods (where applicable);
- Description of any modifications to the test procedure;
- Indication of controls used for direct MTT-reducers and/or colouring test chemicals, if applicable;
- Number of tissue replicates used per test chemical and controls (positive control, negative control, NSMTT, NSCliving and NSCKilled, if applicable);

Results:

- Tabulation of data from individual test chemicals and control substances for each run (including repeat experiments where applicable) and each replicate measurement, including OD value or FD peak area, percent tissue viability, mean percent tissue viability, difference between tissue replicates or SD, and final prediction;
- If applicable, results of controls used for direct MTT-reducers and/or coloured test chemicals, including OD value or FD peak area, %NSMTT, %NSCliving, %NSCKilled, difference between tissue replicates or SD, final correct percent tissue viability, and final prediction;

- Results obtained with the test chemical(s) and control substances in relation to the define run and test acceptance criteria;
- Description of other effects observed, e.g., colouration of the tissues by a coloured test chemical;

Discussion of the results

Conclusion

LITERATURE

- (1) United Nations (UN) (2021), Globally Harmonized System of Classification and Labelling of Chemicals (GHS). Ninth revised edition, New York and Geneva, United Nations Publications. Available at: [<https://unece.org/transport/standards/transport/dangerous-goods/ghs-rev9-2021>].
- (2) OECD (2017). Guideline for Testing of Chemicals No. 405: Acute Eye Irritation/Corrosion. Organisation for Economic Cooperation and Development, Paris. Available at: [<http://www.oecd.org/env/testguidelines>].
- (3) OECD (2018). Guidance Document on an Integrated Approach on Testing and Assessment for Serious Eye Damage and Eye irritation. Series on Testing and Assessment No.263. ENV Publications, Organisation for Economic Cooperation and Development, Paris.
- (4) OECD (2020). Guideline for Testing of Chemicals No. 437: Bovine Corneal Opacity and Permeability Test Method for Identifying i) Chemicals Inducing Serious Eye Damage and ii) Chemicals Not Requiring Classification for Eye Irritation or Serious Eye Damage. Organisation for Economic Cooperation and Development, Paris. Available at: [<http://www.oecd.org/env/testguidelines>].
- (5) OECD (2018). Guideline for Testing of Chemicals No. 438: Isolated Chicken Eye Test Method for Identifying i) Chemicals Inducing Serious Eye Damage and ii) Chemicals Not Requiring Classification. Organisation for Economic Cooperation and Development, Paris. Available at: [<http://www.oecd.org/env/testguidelines>].
- (6) OECD (2017). Guideline for Testing of Chemicals No. 460: Fluorescein Leakage Test Method for Identifying Ocular Corrosives and Severe Irritants. Organisation for Economic Co-operation and Development, Paris. Available at: [<http://www.oecd.org/env/testguidelines>].
- (7) OECD (2020). Guideline for Testing of Chemicals No. 491: Short Time Exposure *In Vitro* Test Method for Identifying i) Chemicals Inducing Serious Eye Damage and ii) Chemicals Not Requiring Classification for Eye Irritation or Serious Eye Damage. Organisation for Economic Co-operation and Development, Paris. Available at: [<http://www.oecd.org/env/testguidelines>].
- (8) OECD (2019). Guideline for Testing of Chemicals No. 492: Short Time Exposure *In Vitro* Test Method for Identifying i) Chemicals Inducing Serious Eye Damage and ii) Chemicals Not Requiring Classification for Eye Irritation or Serious Eye Damage. Organisation for Economic Co-operation and Development, Paris. Available at: [<http://www.oecd.org/env/testguidelines>].
- (9) OECD (2021). Guideline for Testing of Chemicals No. 494: Vitrigel-Eye Irritancy Test Method for Identifying Chemicals Not Requiring Classification and Labelling for Eye Irritation or Serious Eye Damage. Organisation for Economic Co-operation and Development, Paris. Available at: [<http://www.oecd.org/env/testguidelines>].
- (10) OECD (2019). Guideline for Testing of Chemicals No. 496: *In vitro* Macromolecular Test Method for Identifying Chemicals Inducing Serious Eye Damage and Chemicals Not Requiring Classification for Eye Irritation or Serious Eye Damage. Organization for Economic Co-operation and Development, Paris. Available at: [<http://www.oecd.org/env/testguidelines>].
- (11) Alépée, N., Leblanc, V., Grandidier, MH., Teluob, S., Tagliati, V., Adriaens, E., Michaut, V. (2020) Development of the SkinEthic HCE Time-to-Toxicity test method for identifying liquid chemicals not requiring classification and labelling and liquids inducing serious eye damage and eye irritation. *Toxicol. in Vitro* 69, 104960.
- (12) Alépée, N., Leblanc, V., Grandidier, MH., Teluob, S., Viricel, A., Adriaens, E., Michaut, V. (2021). SkinEthic HCE Time-to-Toxicity on Solids: A test method for distinguishing chemicals inducing serious eye damage, eye irritation and not requiring classification and labelling. *Toxicol. in Vitro* 75, 105203.

- (13) Alépée, N., Grandidier, M.H., Teluob, S., Amaral, F., Caviola, E., De Servi, B., Martin, S., Meloni, M., Nardelli, L., Padelou, C., Tagliati, V., Viricel, A., Adriaens, E., Michaut, V. (2022). Validation of the SkinEthic HCE Time-to-Toxicity test method for eye hazard classification of chemicals according to UN GHS. *Toxicol. in Vitro*, 80, 105319. <https://doi.org/10.1016/j.tiv.2022.105319>
- (14) OECD (2022). Peer Review Panel Report on the scientific validity of the SkinEthic™ Human Corneal Epithelium (HCE) Time-to-Toxicity Test (TTT) on liquids (TTL) and on solids (TTS). Series on Testing and Assessment, No. 355. Organisation for Economic Cooperation and Development, Paris. Available at: [<https://www.oecd.org/chemicalsafety/testing/series-testing-assessment-publications-number.htm>]
- (15) Mosmann, T. (1983). Rapid Colorimetric Assay for Cellular Growth and Survival: Application to Proliferation and Cytotoxicity Assays. *J. Immunol. Methods* 65, 55-63.
- (16) Draize, J.H., Woodard, G., Calvery, H.O. (1944). Methods for the Study of Irritation and Toxicity of Substances Applied Topically to the Skin and Mucous Membranes. *Journal of Pharmacol. and Exp. Therapeutics* 82, 377-390
- (17) Nguyen, D.H., Beuerman, R.W., De Wever, B., Rosdy, M. (2003). Three-dimensional construct of the human corneal epithelium for *in vitro* toxicology. In: Salem, H., Katz, S.A. (Eds), *Alternative Toxicological Methods*, CRC Press, pp. 147-159.
- (18) Meloni, M., De Servi, B., Marasco, D., Del Prete, S. (2011). Molecular mechanism of ocular surface damage: Application to an *in vitro* dry eye model on human corneal epithelium. *Molecular Vision* 17, 113-126.
- (19) Scott, L., Eskes, C., Hoffmann, S., Adriaens, E., Alépée, N., Bufo, M., Clothier, R., Facchini, D., Faller, C., Guest, R., Harbell, J., Hartung, T., Kamp, H., Le Varlet, B., Meloni, M., McNamee, P., Osborne, R., Pape, W., Pfannenbecker, U., Prinsen, M., Seaman, C., Spielman, H., Stokes, W., Trouba, K., Van den Berghe, C., Van Goethem, F., Vassallo, M., Vinardell, P., Zuang, V. (2010). A Proposed Eye Irritation Testing Strategy to Reduce and Replace *In Vivo* Studies Using Bottom-Up and Top-Down Approaches. *Toxicol. In Vitro* 24, 1-9.
- (20) Hackett, R.B., McDonald, T.O. (1991). Eye Irritation. In *Advances in Modern Toxicology: Dermatotoxicology* Marzulli F.N. and Maibach H.I. (Eds.), 4th Edition, pp. 749–815. Washington, DC, USA: Hemisphere Publishing Corporation.
- (21) Fox, D.A., Boyes, W.K. (2008). Toxic Responses of the Ocular and Visual System. In *Casarett and Doull's Toxicology: The Basic Science of Poisons* Klaassen C.D.(Ed.), 7th Edition, pp. 665–697. Withby, ON, Canada: McGraw-Hill Ryerson.
- (22) Jester, J.V., Li, H.F., Petroll, W.M., Parker, R.D., Cavanagh, H.D., Carr, G.J., Smith, B., Maurer, J.K. (1998). Area and Depth of Surfactant Induced Corneal Injury Correlates with Cell Death. *Invest. Ophthalmol. Vis. Sci.* 39, 922–936.
- (23) Maurer, J.K., Parker, R.D., Jester, J.V. (2002). Extent of Corneal Injury as the Mechanistic Basis for Ocular Irritation: Key Findings and Recommendations for the Development of Alternative Assays. *Reg. Tox. Pharmacol.* 36, 106-117.
- (24) Jester, J.V., Li, L., Molai, A., Maurer, J.K. (2001). Extent of Corneal Injury as a Mechanistic Basis for Alternative Eye Irritation Tests. *Toxicol. In Vitro* 15, 115–130.
- (25) Jester, J.V., Petroll, W.M., Bean, J., Parker, R.D., Carr, G.J., Cavanagh, H.D., Maurer, J.K. (1998). Area and Depth of Surfactant-Induced Corneal Injury Predicts Extent of Subsequent Ocular Responses. *Invest. Ophthalmol. Vis. Sci.* 39, 2610–2625.
- (26) Adriaens, E., Barroso, J., Eskes, C., Hoffmann, S., McNamee, P., Alépée, N., Bessou-Touya, S., De Smedt, A., De Wever, B., Pfannenbecker, U., Tailhardat, M., Zuang, V. (2014). Retrospective Analysis of the Draize Test for Serious Eye Damage/Eye Irritation: Importance of Understanding the *in vivo* Endpoints Under UN GHS/EU CLP for the Development and Evaluation of *In Vitro* Test

Methods. Arch. Toxicol. 88, 701-723.

- (27) Barroso, J., Pfannenbecker, U., Adriaens, E., Alépée, N., Cluzel, M., De Smedt, A., Hibatallah, J., Klaric, M., Mewes, K.R., Millet, M., Templier, M., McNamee, P. (2017). Cosmetics Europe compilation of historical serious eye damage/eye irritation in vivo data analysed by drivers of classification to support the selection of chemicals for development and evaluation of alternative methods/strategies: the Draize eye test Reference Database (DRD). Arch. Toxicol. 91, 521-547.
- (28) SkinEthic™ HCE Time-to-Toxicity on Liquids SOP, Version 1. (June 16, 2020). *In Vitro* Eye Irritation Time-To-Toxicity Test Method On Liquids (TTL): SkinEthic™ Human Corneal Epithelium Model (HCE). Available at: [www.episkin.com/eye-irritation/Time-To-Toxicity].
- (29) SkinEthic™ HCE Time-to-Toxicity on Solids SOP, Version 1. (June 16, 2020). *In Vitro* Eye Irritation Time-To-Toxicity Test Method On Solids (TTS): SkinEthic™ Human Corneal Epithelium Model (HCE). Available at: [www.episkin.com/eye-irritation/Time-To-Toxicity].
- (30) Alépée, N., Barroso, J., De Smedt, A., De Wever, B., Hibatallah, J., Klaric, M., Mewes, K.R., Millet, M., Pfannenbecker, U., Tailhardat, M., Templier, M., McNamee, P. (2015). Use of HPLC/UPLC-Spectrophotometry for Detection of Formazan in *In Vitro* Reconstructed Human Tissue (RhT)-Based Test Methods Employing the MTT-Reduction Assay to Expand their Applicability to Strongly Coloured Test Chemicals. Toxicol. In Vitro 29, 741-761
- (31) US FDA (2018). Guidance for Industry: Bioanalytical Method Validation. U.S. Department of Health and Human Services, Food and Drug Administration. May 2018. Available at: [<http://www.fda.gov/downloads/Drugs/Guidances/ucm070107.pdf>].
- (32) OECD (2005). Series on Testing and Assessment No. 34: Guidance Document on the Validation and International Acceptance of New or Updated Test Methods for Hazard Assessment. Organisation for Economic Cooperation and Development, Paris. Available at [https://www.oecd.org/officialdocuments/publicdisplaydocumentpdf/?cote=env/jm/mono\(2005\)14&doclanguange=en](https://www.oecd.org/officialdocuments/publicdisplaydocumentpdf/?cote=env/jm/mono(2005)14&doclanguange=en)

ANNEX 1: DEFINITIONS

Balanced Accuracy: Average of the proportion of the correct predictions of each category. The closeness of agreement between test method results and accepted reference values taking into account imbalance of the class distribution.

Benchmark chemical: A chemical used as a standard for comparison to a test chemical. A benchmark chemical should have the following properties: (i) consistent and reliable source(s) for its identification and characterisation; (ii) structural, functional and/or chemical or product class similarity to the chemical(s) being tested; (iii) known physicochemical characteristics; (iv) supporting data on known effects; and (v) known potency in the range of the desired response.

Chemical: A substance or mixture.

Concordance: See "Accuracy".

Cornea: The transparent part of the front of the eyeball that covers the iris and pupil and admits light to the interior.

CV: Coefficient of Variation.

Dev: Deviation.

Eye irritation: Production of changes in the eye, which are fully reversible, occurring after the exposure of the eye to a substance or mixture. Interchangeable with "Reversible effects on the eye" and with "UN GHS Category 2" (1) with the optional Categories 2A (effects fully reversible within 21 days) and 2B (effects fully reversible within 7 days).

Formazan dye (FD): Chromogenic product of the reduction of MTT.

Hazard: Inherent property of an agent or situation having the potential to cause adverse effects when an organism, system or (sub) population is exposed to that agent.

HCE: SkinEthic™ Human Corneal Epithelium

HPLC: High Performance Liquid Chromatography.

IC₅₀: Concentration at which a benchmark chemical reduces the viability of the tissues by 50% following a fixed exposure time (e.g., 30 minutes-treatment with SDS).

Irreversible effects on the eye: See "Serious eye damage".

LLOQ: Lower Limit of Quantification.

LogP: Logarithm of the octanol-water partition coefficient

Mixture: A mixture or a solution composed of two or more substances in which they do not react (1).

Mono-constituent substance: A substance, defined by its quantitative composition, in which one main constituent is present to at least 80% (w/w).

Multi-constituent substance: A substance, defined by its quantitative composition, in which more than one main constituent is present in a concentration $\geq 10\%$ (w/w) and $< 80\%$ (w/w). A multi-constituent substance is the result of a manufacturing process. The difference between mixture and multi-constituent substance is that a mixture is obtained by blending of two or more substances without chemical reaction. A multi-constituent substance is the result of a chemical reaction.

MTT: 3-(4,5-Dimethylthiazol-2-yl)-2,5-diphenyltetrazolium bromide; Thiazolyl blue tetrazolium bromide. (CAS RN 298-93-1)

Negative control: A sample containing all components of a test system and treated with a substance known not to induce a positive response in the test system. This sample is processed with test chemical-treated samples and other control samples and is used to determine 100% tissue viability.

Not Classified: Chemicals that are not classified for Eye irritation (UN GHS Category 2, 2A, or 2B) or Serious eye damage (UN GHS Category 1). Interchangeable with “UN GHS No Category”.

NSC_{killed}: Non-Specific Colour in killed tissues. **NSC_{living}:** Non-Specific Colour in living tissues. **NSMTT:** Non-Specific MTT reduction.

OD: Optical Density.

Positive control: A sample containing all components of a test system and treated with a substance known to induce a positive response in the test system. This sample is processed with test chemical-treated samples and other control samples. To ensure that variability in the positive control response across time can be assessed, the magnitude of the positive response should not be excessive.

Relevance: Description of relationship of the test to the effect of interest and whether it is meaningful and useful for a particular purpose. It is the extent to which the test correctly measures or predicts the biological effect of interest. Relevance incorporates consideration of the accuracy (concordance) of a test method (32).

Reliability: Measures of the extent that a test method can be performed reproducibly within and between laboratories over time, when performed using the same protocol. It is assessed by calculating intra- and inter-laboratory reproducibility and intra-laboratory repeatability (32).

Replacement test: A test which is designed to substitute for a test that is in routine use and accepted for hazard identification and/or risk assessment, and which has been determined to provide equivalent or improved protection of human or animal health or the environment, as applicable, compared to the accepted test, for all possible testing situations and chemicals (32).

Reproducibility: The agreement among results obtained from repeated testing of the same test chemical using the same test protocol (See "Reliability") (32).

Reversible effects on the eye: See “Eye irritation”.

RhCE: Reconstructed human Cornea-like Epithelium.

Run: A run consists of one or more test chemicals tested concurrently with a negative control and with a positive control.

SD: Standard Deviation.

Serious eye damage: refers to the production of tissue damage in the eye, or serious physical decay of vision, which is not fully reversible, occurring after exposure of the eye to a substance or mixture. Interchangeable with “Irreversible effects on the eye” and with “UN GHS Category 1” (1).

Standard Operating Procedures (SOP): Formal, written procedures that describe in detail how specific routine, and test-specific, laboratory operations should be performed. They are required by GLP.

Substance: Chemical elements and their compounds in the natural state or obtained by any production process, including any additive necessary to preserve the stability of the product and any impurities deriving from the process used, but excluding any solvent which may be separated without affecting the stability of the substance or changing its composition (1).

Test: A single test chemical concurrently tested in a minimum of two tissue replicates as defined in the corresponding SOP.

Tissue viability: Parameter measuring total activity of a cell population in a reconstructed tissue as their ability to reduce the vital dye MTT, which, depending on the endpoint measured and the test design used, correlates with the total number and/or vitality of living cells.

Test chemical: The term "test chemical" is used to refer to what is being tested.

TTL: SkinEthic™ HCE Time-to-Toxicity for Liquids chemicals

TTS: SkinEthic™ HCE Time-to-Toxicity for Solids chemicals

TTT: SkinEthic™ HCE Time-to-Toxicity test method

ULOQ: Upper Limit of Quantification.

United Nations Globally Harmonized System of Classification and Labelling of Chemicals (UN GHS): A system proposing the classification of chemicals (substances and mixtures) according to standardised types and levels of physical, health and environmental hazards, and addressing corresponding communication elements, such as pictograms, signal words, hazard statements, precautionary statements and safety data sheets, so that to convey information on their adverse effects with a view to protect people (including employers, workers, transporters, consumers and emergency responders) and the environment (1).

UN GHS Category 1: See "Serious eye damage".

UN GHS Category 2: See "Eye irritation".

UN GHS No Category: Chemicals that do not meet the requirements for classification as UN GHS Category 1 or 2 (2A or 2B). Interchangeable with "Not Classified".

UPLC: Ultra-High Performance Liquid Chromatography.

UVCB: substances of unknown or variable composition, complex reaction products or biological materials.

Valid test method: A test method considered to have sufficient relevance and reliability for a specific purpose and which is based on scientifically sound principles. A test method is never valid in an absolute sense, but only in relation to a defined purpose (31).

Validated test method: A test method for which validation studies have been completed to determine the relevance (including accuracy) and reliability for a specific purpose. It is important to note that a validated test method may not have sufficient performance in terms of accuracy and reliability to be found acceptable for the proposed purpose (31).

VRM: Validated Reference Method.

Weight-of-evidence: The process of considering the strengths and weaknesses of various pieces of information in reaching and supporting a conclusion concerning the hazard potential of a test chemical.

ANNEX 2: MAIN TEST METHOD COMPONENTS OF THE SKINETHIC™ HCE TTT TEST METHOD VALIDATED FOR EYE HAZARD IDENTIFICATION

Test Method Component		SkinEthic™ HCE TTT (VRM)	
Protocol¹	TTL protocol: Liquid and viscous test chemicals (pipetteable)	TTS protocol: Solid test chemicals (not pipetteable)	
Model surface	0.5 cm ²	0.5 cm ²	
Number of tissue replicates	At least 2	At least 2	
Pre- or post-check for colour interference	10 µL + 90 µL H ₂ O mixed for 30±2 min at Room Temperature (RT) ² → If test chemical is coloured, living adapted controls should be performed	10 mg + 90 µL H ₂ O mixed for 30±2 min at RT → If test chemical is coloured, living adapted controls should be performed	
Pre- or post-check for direct tetrazolium reduction	50-80 µL chemical + 300 µL MTT 1 mg/mL solution for 180±15 min at 37±2°C, 5±1% CO ₂ , ≥95% RH →if solution turns blue/purple, water-killed adapted controls should be performed (50-80 µL of sterile deionized water in MTT solution is used as negative control)	30-80 mg chemical + 300 µL MTT 1 mg/mL solution for 180±15 min at 37±2°C, 5±1% CO ₂ , ≥95% RH →if solution turns blue/purple, water-killed adapted controls should be performed (30-80 µL of sterile deionized water in MTT solution is used as negative control)	
Treatment doses and application	10 µL Ca ²⁺ /Mg ²⁺ -free DPBS + 80±2 µL (160 µL/cm ²)	80±2 mg (160 mg/cm ²) (grounded if needed to cover the whole tissue surface)	

¹ For details of the protocol, the latest version of the Standard Operating Procedures (SOP) should be consulted and employed when implementing the test method in a laboratory (28-29).

² Room Temperature (RT): 18-28°C

OECD/OCDE 492B


		For viscous, use a nylon mesh	+ 80 μ L distilled water Use a nylon mesh
Exposure time and temperature		5 min (\pm 15 sec) (neat), 16 min (\pm 1 min) (at 20 % (w/v) in distilled water) in culture medium at RT and 120 min (\pm 2 min) (at 20 % (w/v) in distilled water) in culture medium at 37 \pm 2 $^{\circ}$ C, 5 \pm 1% CO ₂ , \geq 95% RH	30 min (\pm 2 min) and 120 min (\pm 5 min) (neat) in culture medium at 37 \pm 2 $^{\circ}$ C, 5 \pm 1% CO ₂ , \geq 95% RH
Rinsing at room temperature		25 mL of Ca ²⁺ /Mg ²⁺ -free DPBS (2 mL per push)	25 mL of Ca ²⁺ /Mg ²⁺ -free DPBS (2 mL per push)
Post-soak immersion		10 min (\pm 1 min) at RT in 4 mL culture medium	30 min (\pm 2 min) at RT in 4 mL culture medium
Negative control		80 \pm 2 μ L Ca ²⁺ /Mg ²⁺ -free DPBS tested concurrently	80 \pm 2 μ L Ca ²⁺ /Mg ²⁺ -free DPBS tested concurrently
Positive control		80 \pm 2 μ L Methyl acetate tested concurrently (CAS RN 79-20-9)	80 \pm 2 μ L Lactic acid 1% in water (w/v) tested concurrently (CAS RN 50-21-5)
Tetrasolium salt solution		300 μ L 1 mg/mL	300 μ L 1 mg/mL
Tetrazolium salt incubation and temperature		180 min (\pm 15 min) at 37 \pm 2 $^{\circ}$ C, 5 \pm 1% CO ₂ , \geq 95% RH	180 min (\pm 15 min) at 37 \pm 2 $^{\circ}$ C, 5 \pm 1% CO ₂ , \geq 95% RH
Extraction solvent		1.5 mL isopropanol (750 μ L under and 750 μ L over) (extraction from top and bottom of insert)	1.5 mL isopropanol (extraction from bottom of insert)
Extraction time and temperature		At least 2 h with shaking (\sim 120 rpm) at RT or at least overnight without shaking at 4-10 $^{\circ}$ C	At least 2 h with shaking (\sim 120 rpm) at RT
OD reading		570 nm (540 - 600 nm) without reference filter	570 nm (540 - 600 nm) without reference filter
Tissue Quality Control		30 \pm 2 min treatment with SDS (50 μ L) 1.0 mg/mL \leq IC ₅₀ \leq 3.2 mg/mL	30 \pm 2 min treatment with SDS (50 μ L) 1.0 mg/mL \leq IC ₅₀ \leq 3.2 mg/mL
Acceptance Criteria		1. Mean OD of the tissue replicates treated with the negative control should be > 1.0 and	1. Mean OD of the tissue replicates treated with the negative control should be > 1.0 and

OECD/OCDE

492B

	<p>≤ 2.5</p> <p>2. Mean viability of the PC, expressed as % of the NC, is ≤ 50% at the time exposure of 5 minutes and > 50% at both times 16 and 120 minutes.</p> <p>3. The difference of viability between two tissue replicates should be less than 20%.</p>	<p>≤ 2.5</p> <p>2. Mean viability of the PC, expressed as % of the NC, is > 40% at the time exposure of 30 minutes and 20% < PC ≤ 70% at the time exposure of 120 minutes.</p> <p>3. The difference of viability between two tissue replicates should be less than 20%.</p>
--	---	---

ANNEX 3: ILLUSTRATIVE FLOWCHART PROVIDING GUIDANCE ON HOW TO IDENTIFY AND HANDLE DIRECT MTT-REDUCERS AND/OR COLOUR INTERFERING CHEMICALS, BASED ON THE SKINETHIC™ HCE TTT SOP



ANNEX IV: KEY PARAMETERS AND ACCEPTANCE CRITERIA FOR QUALIFICATION OF AN HPLC/UPLC-SPECTROPHOTOMETRY SYSTEM FOR MEASUREMENT OF MTT FORMAZAN EXTRACTED FROM RhCE TISSUE CONSTRUCTS

Parameter	Protocol Derived from FDA Guidance (43)(45)	Acceptance Criteria
Selectivity	Analysis of isopropanol, living blank (isopropanol extract from living RhCE tissue constructs without any treatment), dead blank (isopropanol extract from killed RhCE tissue constructs without any treatment), and of a dye (e.g., methylene blue)	$\text{Area}_{\text{interference}} \leq 20\% \text{ of Area}_{\text{LLOQ}}^1$
Precision	Quality Controls (i.e., MTT formazan at 1.6 $\mu\text{g/mL}$, 16 $\mu\text{g/mL}$ and 160 $\mu\text{g/mL}$) in isopropanol (n=5)	$\text{CV} \leq 15\% \text{ or } \leq 20\% \text{ for the LLOQ}$
Accuracy	Quality Controls in isopropanol (n=5)	$\% \text{Dev} \leq 15\% \text{ or } \leq 20\% \text{ for LLOQ}$
Matrix Effect	Quality Controls in living blank (n=5)	$85\% \leq \% \text{Matrix Effect} \leq 115\%$
Carryover	Analysis of isopropanol after an ULOQ ² standard	$\text{Area}_{\text{interference}} \leq 20\% \text{ of Area}_{\text{LLOQ}}$
Reproducibility (intra-day)	3 independent calibration curves (based on 6 consecutive 1/3 dilutions of MTT formazan in isopropanol starting at ULOQ, i.e., 200 $\mu\text{g/mL}$); Quality Controls in isopropanol (n=5)	Calibration Curves: $\% \text{Dev} \leq 15\% \text{ or } \leq 20\% \text{ for LLOQ}$
Reproducibility (inter-day)	Day 1: 1 calibration curve and Quality Controls in isopropanol (n=3) Day 2: 1 calibration curve and Quality Controls in isopropanol (n=3) Day 3: 1 calibration curve and Quality Controls in isopropanol (n=3)	Quality Controls: $\% \text{Dev} \leq 15\% \text{ and } \text{CV} \leq 15\%$
Short Term Stability of MTT Formazan in RhCE Tissue Extract	Quality Controls in living blank (n=3) analysed the day of the preparation and after 24 hours of storage at room temperature	$\% \text{Dev} \leq 15\%$
Long Term Stability of MTT Formazan in RhCE Tissue Extract, if required	Quality Controls in living blank (n=3) analysed the day of the preparation and after several days of storage at -20°C	$\% \text{Dev} \leq 15\%$

¹LLOQ: Lower Limit of Quantification, defined to cover 1-2% tissue viability, i.e., 0.8 $\mu\text{g/mL}$.

²ULOQ: Upper Limit of Quantification, defined to be at least two times higher than the highest expected MTT formazan concentration in isopropanol extracts from negative controls ($\sim 70 \mu\text{g/mL}$ in the VRM), i.e., 200 $\mu\text{g/mL}$.

ANNEX V: PERFORMANCE OF THE SKINETHIC™ HCE TTT FOR THE LIQUIDS AND SOLIDS USING THE 3X3 MATRIX SHOWING CORRECT, UNDER- AND OVER- PREDICTIONS PER UN GHS CATEGORY

Table 1: Performance of the SkinEthic™ HCE TTL protocol

47.

UN GHS categories	SkinEthic™ HCE TTL - Predicted categories		
	Cat 1	Cat 2	No Cat
Cat 1 (N=21)	85.4%	14.6%	0.0%
Cat 2 (N=25)	20.2%	79.8%	0.0%
No Cat (N=24)	0.0%	20.8%	79.2%

Table 2: Performance of the SkinEthic™ HCE TTS protocol

48.

UN GHS categories	SkinEthic™ HCE TTS - Predicted categories		
	Cat 1	Cat 2	No Cat
Cat 1 (N=29)	74.7%	25.3%	0.0%
Cat 2 (N=19)	15.8%	55.3%	28.9%
No Cat (N=33)	3.0%	25.3%	71.7%

