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ADRA: Amino acid Derivative Reactivity Assay

AOP: Adverse Outcome Pathway

ARE: Antioxidant Responsive Element

ATP: Adenosine Triphosphate

BA: Balanced Accuracy

BR: Borderline Range

BrdU: Bromodeoxyuridine

CD: Cluster of Differentiation

DA: Defined Approach

DASS: Defined Approach for Skin Sensitisation

DIP: Data Interpretation Procedure

DPRA: Direct Peptide Reactivity Assay

ECHA: European Chemicals Agency

EU: European Union

GL: Guideline

h-CLAT: human Cell Line Activation Test

ITS: Integrated Testing Strategy

KE: Key Event

LLNA: Local Lymph Node Assay

MIT: Minimum Induction Threshold

Nrf2: nuclear factor-erythroid 2-related factor 2

NC: Not classified

OECD: Organisation for Economic Co-operation and Development

QSAR: Quantitative Structure-Activity Relationship

QMRF: QSAR Model Reporting Format

REACH: Registration, Evaluation, Authorisation and Restriction of Chemicals

RFI: Relative Fluorescence Intensity

RI: Radio Isotope

TG: Test Guideline

UN GHS: United Nations Globally Harmonized System of Classification and Labelling of Chemicals

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DEREK: Derek Nexus

OECD TB: OECD QSAR Toolbox

Note: UN GHS 1A correspond to strong sensitizers and UN GHS 1B correspond to other (moderate to weak) sensitizers. Not classified are considered non-sensitizers. *Cysteine-only depletion thresholds are used in the case of co-elution with the lysine peptide.

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DA Performance vs. LLNA Data (N=134)	2o3
Accuracy (%)	83%
Sensitivity (%)	82%
Specificity (%)	85%
Balanced Accuracy (%)	84%

Note: Accuracy is the correct classification rate, sensitivity is the true positive rate, specificity is the true negative rate, and balanced accuracy is the average of sensitivity and specificity. Performance is reported based on DPRA, KeratinoSens™, and h-CLAT. Statistics reflect conclusive predictions only; inconclusive predictions are shown in grey. Additional performance characterisation is available in the Supporting document to the Guideline (GL) on Defined Approaches (DAs) for Skin Sensitisation.

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DA Performance vs. Human Data (N=55)	2o3
Accuracy (%)	89%
Sensitivity (%)	89%
Specificity (%)	88%
Balanced Accuracy (%)	88%

Note: Accuracy is the correct classification rate, sensitivity is the true positive rate, specificity is the true negative rate, and balanced accuracy is the average of sensitivity and specificity with respect to HPPT data. Performance is reported based on DPRA, KeratinoSens™, and h-CLAT. Statistics reflect conclusive predictions only; inconclusive predictions are shown in grey. Additional performance characterisation is available in the Supporting document to the Guideline (GL) on Defined Approaches (DAs) for Skin Sensitisation (1).

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DA Performance vs. LLNA Data (N=159)	ITSv1
Accuracy (%)	87%
Sensitivity (%)	92%
Specificity (%)	70%
Balanced Accuracy (%)	81%

Note: Accuracy is the correct classification rate, sensitivity is the true positive rate, specificity is the true negative rate, and balanced accuracy is the average of sensitivity and specificity with respect to LLNA data. Statistics reflect high confidence predictions only; inconclusive predictions are shown in grey. Additional performance characterisation is available in the Supporting document to the Guideline (GL) on Defined Approaches (DAs) for Skin Sensitisation (1).

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DA Performance vs. Human Data (N=64)	ITSv1
Accuracy (%)	86%
Sensitivity (%)	93%
Specificity (%)	44%
Balanced Accuracy (%)	69%

Note: Accuracy is the correct classification rate, sensitivity is the true positive rate, specificity is the true negative rate, and balanced accuracy is the average of sensitivity and specificity with respect to Human HPPT-based data. Statistics reflect conclusive predictions only; inconclusive predictions are shown in grey. Additional performance characterisation is available in the Supporting document to the Guideline (GL) on Defined Approaches (DAs) for Skin Sensitisation (1).

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DA Performance vs. LLNA Data (N=156)	ITSv2
Accuracy (%)	88%
Sensitivity (%)	93%
Specificity (%)	67%
Balanced Accuracy (%)	80%

Note: Accuracy is the correct classification rate, sensitivity is the true positive rate, specificity is the true negative rate, and balanced accuracy is the average of sensitivity and specificity with respect to LLNA data. Statistics reflect conclusive predictions only; inconclusive predictions are shown in grey. Additional performance characterisation is available in the Supporting document to the Guideline (GL) on Defined Approaches (DAs) for Skin Sensitisation (1).

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DA Performance vs. Human Data (N=62)	ITSv2
Accuracy (%)	87%
Sensitivity (%)	94%
Specificity (%)	44%
Balanced Accuracy (%)	69%

Note: Accuracy is the correct classification rate, sensitivity is the true positive rate, specificity is the true negative rate, and balanced accuracy is the average of sensitivity and specificity with respect to Human HPPT-based data. Statistics reflect conclusive predictions only; inconclusive predictions are shown in grey. Additional performance characterisation is available in the Supporting document to the Guideline (GL) on Defined Approaches (DAs) for Skin Sensitisation(1).

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DA/Method	Information Sources	Capability (Hazard and/or Potency)	Hazard Performance vs. LLNA	Hazard Performance vs. Human	Potency Performance vs. LLNA (Accuracy)	Potency Performance vs. Human (Accuracy)
2o3 DA	DPRA, KeratinoSens™, h-CLAT	Hazard	84% BA, 82% Sens, 85% Spec	88% BA, 89% Sens, 88% Spec	-	-
ITSv1 DA	DPRA, h-CLAT, DEREK Nexus v6.1.0	Hazard, Potency	81% BA, 92% Sens, 70% Spec	69% BA, 93% Sens, 44% Spec	70% NC, 71% 1B, 74% 1A	44% NC, 77% 1B, 65% 1A
ITSv2 DA	DPRA, h-CLAT, OECD QSAR Toolbox v4.5	Hazard, Potency	80% BA, 93% Sens, 67% Spec	69% BA, 94% Sens, 44% Spec	67% NC, 72% 1B, 72% 1A	44% NC, 80% 1B, 67% 1A
LLNA (provided for comparison)	<i>in vivo</i>	Hazard, Potency	-	58% BA, 94% Sens, 22% Spec	-	25% NC, 74% 1B, 56% 1A

Note: For hazard performance, sensitivity (Sens) is the true positive rate, specificity (Spec) is the true negative rate, and balanced accuracy (BA) is the average of sensitivity and specificity. Due to the imbalanced nature of the reference data, the measures of specificity are more uncertain than the measures of sensitivity. For potency performance, accuracy reflects correct classification rate within each UN GHS sub-category. Due to the imbalanced nature of the reference data, the measures of accuracy are more uncertain for smaller classes, e.g. for NC chemicals. Statistics reflect conclusive DA predictions only. This represents the data available at the time of initial guideline adoption.

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^an: number of test chemicals out of the 24 chemicals assessed in the validation study for which at least three test runs were available.

^bPublished in Gabbert *et al.* (2020).

For the three laboratories participating in the validation study and only considering test chemicals for which at least three test runs were available (i.e. 13 to 14), the borderline range around the 6.38% mean peptide depletion cut-off varied between 4.54 to 5.49 (lower boundary) and 7.42 to 9.08% (upper boundary). *The DPRA mean borderline range of all participating laboratories was 4.95 to 8.32%*, which was comparable to the range of 5.29 to 7.69% derived from historical data of a routine testing lab (assessing 385 substances) and to the experimentally determined range of 5.45 to 7.31% when repeatedly testing a single substance (i.e. EGDMA) (the latter two published in Gabbert *et al.*, 2020). Likewise, the different borderline ranges determined for the cysteine-only depletion model were also very comparable with each other.

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^a n: number of test chemicals out of the 28 chemicals assessed in the validation study for which at least three test runs were available.

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Table 2.4. h-CLAT borderline ranges determined based on the log pooled median absolute deviations.

Data source	RFI CD54 (cut-off 200)	RFI CD86 (cut-off 150)
Validation study lab 1 (n ^a = 24)	152 - 264	125 - 181
Validation study lab 2 (n ^a = 24)	153 - 261	125 - 181
Validation study lab 3 (n ^a = 24)	161 - 248	115 - 196
Validation study lab 4 (n ^a = 24)	162 - 247	125 - 180
Validation study mean (lab 1-4)	157 - 255	122 - 184
BASF SE historical data ^b (n = 136)	170 - 235	132 - 170

^a n: number of test chemicals of the 24 chemicals assessed in the ring trial for which at least

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- 1) OECD (2022). OECD Guidelines for the Testing of Chemicals No.406. Skin Sensitisation Guinea Pig Maximisation Test and Buehler Test, Organisation for Economic Cooperation and Development, Paris. Available at: ([oecd-ilibrary.org](https://www.oecd-ilibrary.org))
- 2) OECD (2010). OECD Guidelines for the Testing of Chemicals No. 429. Skin Sensitization: Local Lymph Node Assay (LLNA), Organisation for Economic Cooperation and Development, Paris. Available at: ([oecd-ilibrary.org](https://www.oecd-ilibrary.org))
- 3) OECD (2010). OECD Guidelines for the Testing of Chemicals No. 442A. Skin Sensitization: Local Lymph Node Assay: DA, Organisation for Economic Cooperation and Development, Paris. Available at: ([oecd-ilibrary.org](https://www.oecd-ilibrary.org))
- 4) OECD (2018). OECD Guidelines for the Testing of Chemicals No. 442B. Local lymph node assay: BrdU-ELISA or -FCM, Organisation for Economic Cooperation and Development, Paris. Available at: ([oecd-ilibrary.org](https://www.oecd-ilibrary.org))
- 5) OECD (2024). OECD Guidelines for the Testing of Chemicals No. 442C: *In chemico* Skin Sensitisation Assays addressing the Adverse Outcome Pathway key event on covalent binding to proteins, Organisation for Economic Cooperation and Development, Paris. Available at: ([oecd-ilibrary.org](https://www.oecd-ilibrary.org))
- 6) OECD (2024). OECD Key event based test Guideline No. 442D: *In vitro* Skin Sensitisation Assays Addressing AOP Key Event on Keratinocyte Activation. Organisation for Economic Cooperation and Development, Paris. Available at: ([oecd-ilibrary.org](https://www.oecd-ilibrary.org))
- 7) OECD (2024). OECD Key event based test Guideline No. 442E: *In vitro* Skin Sensitisation Assays Addressing the Key Event on Activation of Dendritic Cells on the Adverse Outcome Pathway for Skin Sensitisation. Organisation for Economic Cooperation and Development, Paris. Available at: ([oecd-ilibrary.org](https://www.oecd-ilibrary.org))

- 8) OECD (2012). Series on Testing and Assessment No. 168. The Adverse Outcome Pathway for Skin Sensitisation Initiated by Covalent Binding to Proteins. Part 1: Scientific Evidence. Organisation for Economic Cooperation and Development, Paris. Available at: (oecd-ilibrary.org)
- 9) OECD (2023). OECD Guidelines No.497. Defined Approaches on Skin Sensitisation, Organisation for Economic Cooperation and Development, Paris. Available at: (oecd-ilibrary.org)
- 10) OECD (2021). Series on Testing and Assessment No. 336. Supporting Document to the OECD Guideline 497 on Defined Approaches for Skin Sensitisation, Organisation for Economic Cooperation and Development, Paris. Available at: (oecd-ilibrary.org)
- 11) Kolle, S.N., et al. (2021). Assessing Experimental Uncertainty in Defined Approaches: Borderline Ranges for *In Chemico* and *In Vitro* Skin Sensitization Methods Determined from Ring trial Data, *Applied in Vitro Toxicology*, 7(3), 102-101. DOI:10.1089/aivt.2121.0003
- 12) ECHA (2008). see “CHAPTER R.6 – QSARS AND GROUPING OF CHEMICALS” in Guidance on Information Requirements and Chemical Safety Assessment. European Chemicals Agency [[Guidance on Information Requirements and Chemical Safety Assessment - ECHA \(europa.eu\)](https://echa.europa.eu)].
- 13) Chilton, M. L., Macmillan, D. S., Steger-Hartmann, T., Hillegass, J., Bellion, P., Vuorinen, A., Etter, S., Smith, B. P. C., White, A., Sterchele, P., De Smedt, A., Glogovac, M., Glowienke, S., O’Brien, D., & Parakhia, R. (2018). Making reliable negative predictions of human skin sensitisation using an *in silico* fragmentation approach. *Regulatory Toxicology and Pharmacology*, 95, 227–235. <https://doi.org/10.1016/j.yrtph.2018.03.015>

Appendix 1 3 of 3, 2 of 2 および 2 of 3 DA の予測性

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	LLNA	
	Non	Sens
3o3		
Non	10	1
Sens	9	80
Inconclusive	14	53

	LLNA	
	Non	Sens
DA 2o3		
Non	22	19
Sens	4	89
Inconclusive	7	27

	3o3
(N=100)	
Accuracy (%)	90
Sensitivity (%)	99
Specificity (%)	53
Balanced Accuracy (%)	76

	2o3 DA
(N=134)	
Accuracy (%)	83
Sensitivity (%)	82
Specificity (%)	85
Balanced Accuracy (%)	84

2o2(1)	LLNA	
DPRAs and KeratinoSens	Non	Sens
Non	19	40
Sens	10	78
Inconclusive	4	17

2o2(2)	LLNA	
h-CLAT and KeratinoSens	Non	Sens
Non	11	26
Sens	9	59
Inconclusive	13	49

DA performance vs. LLNA data	DPRAs and KeratinoSens
(N=147)	
Accuracy (%)	66
Sensitivity (%)	66
Specificity (%)	66
Balanced Accuracy (%)	66

DA performance vs. LLNA data	h-CLAT and KeratinoSens
(N=105)	
Accuracy (%)	67
Sensitivity (%)	69
Specificity (%)	55
Balanced Accuracy (%)	62

2o2(3)	LLNA	
h-CLAT and DPRAs	Non	Sens
NC	12	24
Sens	9	62
Inconclusive	12	48

DA performance vs. LLNA data	DPRAs and h-CLAT
(N=107)	
Accuracy (%)	69
Sensitivity (%)	72
Specificity (%)	57
Balanced Accuracy (%)	65

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DA performance vs. Human data	DPRAs and KeratinoSens
(N=54)	
Accuracy (%)	74
Sensitivity (%)	74
Specificity (%)	75
Balanced Accuracy (%)	74

DA performance vs. Human data	h-CLAT and KeratinoSens
(N=45)	
Accuracy (%)	67
Sensitivity (%)	69
Specificity (%)	50
Balanced Accuracy (%)	60

2o2(3)	Human	
h-CLAT, DPRAs	NC	Sens
NC	2	10
Sens	4	31
Borderline negative, Borderline positive, Inconclusive, NA	5	12

NA: Not Available

DA performance vs. Human data	DPRAs and h-CLAT
(N=47)	
Accuracy (%)	70
Sensitivity (%)	76
Specificity (%)	33
Balanced Accuracy (%)	54

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Note: Chemicals highlighted in light green have false negatives and those highlighted in yellow have false positives.

LLNA.GHS.BIN: LLNA Binary hazard reference classification

LLNA.GHS.SUB: LLNA Potency reference subcategorisation

2 of 3 DA Call.Conf: 2 of 3 DA Hazard prediction considering confidence workflow

HU.GHS.BIN: Human Binary hazard reference classification

Appendix 3 ITS DA の予測性

1)ITSv1 DA

ITSv1 DA in comparison to LLNA
(GL497 と同じ)

ITSv1 DA	LLNA		
	NC	1B	1A
NC	21	11	0
1B	9	55	10
1A	0	12	28
Inconclusive	3	7	0

71% correct classification overall for potency
80% balanced accuracy overall for hazard
(Sensitivity 91% Specificity 70%)

ITSv1 DA in comparison to Human
(GL497 と同じ)

ITSv1 DA	Human		
	NC	1B	1A
NC	4	4	0
1B	5	24	7
1A	0	3	13
Inconclusive	2	0	1

68% correct classification overall for potency
68% balanced accuracy overall for hazard
(Sensitivity 92% Specificity 44%)

2) ITSv2 DA

ITSv2 DA in comparison to LLNA
(GL497 と同じ)

ITSv2 DA	LLNA		
	NC	1B	1A
NC	20	9	0
1B	10	54	10
1A	0	12	26
Inconclusive	3	10	2

71% correct classification overall for potency
79% balanced accuracy overall for hazard
(Sensitivity 92% Specificity 67%)

ITSv1 DA Score h-CLAT and DPRA in
comparison to LLNA (*in silico* を含まず、
Borderline 未使用)

h-CLAT, DPRA	LLNA		
	NC	1B	1A
NC	22	17	0
1B	8	55	20
1A	0	6	18

65% correct classification overall for potency
79% balanced accuracy overall for hazard
(Sensitivity 85% Specificity 73%)

ITSv1 DA Score h-CLAT and DPRA in
comparison to Human (*in silico* を含まず、
Borderline 未使用)

h-CLAT, DPRA	Human		
	NC	1B	1A
NC	4	7	2
1B	5	22	12
1A	0	1	6

54% correct classification overall for potency
63% balanced accuracy overall for hazard
(Sensitivity 82% Specificity 44%)

ITSv2 DA Score h-CLAT and DPRA in
comparison to LLNA (*in silico* を含まず、
Borderline 未使用)

h-CLAT, DPRA	LLNA (n=141)		
	NC	1B	1A
NC	22	14	0
1B	8	55	18
1A	0	6	18

67% correct classification overall for potency
80% balanced accuracy overall for hazard
(Sensitivity 87% Specificity 73%)

ITSv2 DA in comparison to Human
(GL497 と同じ)

ITSv2 DA	Human		
	NC	1B	1A
NC	4	3	0
1B	5	24	6
1A	0	3	12
Inconclusive	2	1	3

70% correct classification overall for potency
69% balanced accuracy overall for hazard
(Sensitivity 94% Specificity 44%)

ITSv2 DA Score h-CLAT and DPRA in
comparison to Human (*in silico* を含まず、
Borderline 未使用)

h-CLAT, DPRA	Human		
	NC	1B	1A
NC	4	6	1
1B	5	22	11
1A	0	1	6

57% correct classification overall for potency
65% balanced accuracy overall for hazard
(Sensitivity 85% Specificity 44%)

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Note: Note: Chemicals highlighted in light green have misclassification and those highlighted in yellow have false positives.

LLNA.GHS.SUB: LLNA Potency reference subcategorization

ITSv1 DA Score: ITSv1 DA total score

ITSv1.DA Pot.Conf: ITSv1 DA Hazard prediction considering confidence workflow

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Note: Note: Chemicals highlighted in light green have misclassification and those highlighted in yellow have false positives.

LLNA.GHS.SUB: LLNA Potency reference subcategorization

ITSv2 DA Score: ITSv2 DA total score

ITSv2.DA Pot.Conf: ITSv2 DA Hazard prediction considering confidence workflow

HU.GHS.BIN: Human Binary hazard reference classification